

Product datasheet for **TP320433M**

CORO6 (NM_032854) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens coronin 6 (CORO6), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220433 representing NM_032854 Red =Cloning site Green =Tags(s) MSRRVVRQSKFRHVFGQAAKADQAYEDIRVSKVTWDSSFCVNPFLAIIVEAGGGGAFIVLPLAKTGRV DKNYPLVTGHTAPVLDIDWCPHNDNVIASASDDTTIMVWQIPDYTPMRNITEPIITLEGHSKRVGILSWH PTARNVLLSAGGDNVWVGTGEVLLSLDDMHPDVIHVCWNSNGSLLATTCKDKTLRIIDPRKGQVV AERFAAHEGMRPMAVFTRQGHIFTTGFRMSQRELGLWDPNNFEEPVALQEMDTSNGVLLPFYDPDSS I VYLCGKGDSSIRYFEITDEPPFVHYLNTFSSKEPQRGMGFMPKRGLDVSKCEIARFYKLHERKCEPIIMT VPRKSDLFQDDLYPDTGPEPALEADEWLSGQDAEPVLISLRDGYVPPKHRELVRTKRNILDVRPPSGPR RSQSASDAPLSQQHTLETLLEEIKALRERVQAQEQRITALENMLCELVDGTD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	52.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: NP_116243
 Locus ID: 84940
 UniProt ID: Q6QEF8
 RefSeq Size: 2450
 Cytogenetics: 17q11.2
 RefSeq ORF: 1416

Product images:


Coomassie blue staining of purified CORO6 protein (Cat# [TP320433]). The protein was produced from HEK293T cells transfected with CORO6 cDNA clone (Cat# [RC220433]) using MegaTran 2.0 (Cat# [TT210002]).