

Product datasheet for **TP320441L**

ENSA (NM_207042) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human endosulfine alpha (ENSA), transcript variant 1, 1 mg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC220441 representing NM_207042
Red=Cloning site **Green**=Tags(s)

MSQKQEEENPAEETGEEKQDTQEKEGILPERAEEAKLKAKYPSLGQKPGGSDFLMKRLQKGDYKSLHWSV
LLCADEMQKYFDSGDYDMAKAKMKNKQLPSAGPDKNLVTDHIPTQDLPQRKSSLVTSKLAGGQVE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 15.1 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_996925](#)

Locus ID: 2029

UniProt ID: [O43768](#)

RefSeq Size: 1300

Cytogenetics: 1q21.3



[View online »](#)


RefSeq ORF: 411

Synonyms: ARPP-19e

Summary: The protein encoded by this gene belongs to a highly conserved cAMP-regulated phosphoprotein (ARPP) family. This protein was identified as an endogenous ligand for the sulfonylurea receptor, ABCC8/SUR1. ABCC8 is the regulatory subunit of the ATP-sensitive potassium (KATP) channel, which is located on the plasma membrane of pancreatic beta cells and plays a key role in the control of insulin release from pancreatic beta cells. This protein is thought to be an endogenous regulator of KATP channels. In vitro studies have demonstrated that this protein modulates insulin secretion through the interaction with KATP channel, and this gene has been proposed as a candidate gene for type 2 diabetes. At least eight alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified ENSA protein (Cat# [TP320441]). The protein was produced from HEK293T cells transfected with ENSA cDNA clone (Cat# [RC220441]) using MegaTran 2.0 (Cat# [TT210002]).