

Product datasheet for TP320503

JAK2 (NM_004972) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human Janus kinase 2 (a protein tyrosine kinase) (JAK2), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220503 representing NM_004972 Red =Cloning site Green =Tags(s)

MGMACLTMTMEGTSTSSIQNGDISGNANSMKQIDPVLQVYLYHSLGKSEADYLTFFPSGEYVAEEICIA
ASKACGIPVYHNMFALMSETERIWYPPNHVHFHIDESTRHNVLYRIRFYFPRWYCSGSNRAYRHGISRGA
EAPLLDDFVMSYLFAQWRHDFVHGWIKVPVTHETQEECLGMAVLDMMRIAKENDQTPLAIYNSISYKTF
PKCIRAKIQDYHILTRKIRYRFRRIQQFSQCKATARNLKLKYLINLETLSAFYTEKFEVKEPGSGPS
GEEIFATIIITGNGGIQWSRGKHKESETLTEQDLQLYCDFPNIIDVSIKQANQEGSNESRVVTIHKQDGK
NLEIELSSLREALSFVSLIDGYRRLTADAHHYLCKEVAPPVLENIQSNCHGPISMDFAISLKKAGNQT
GLYVLRCSPKDFNKYFLTFAVERENVIEYKHCLITKNENEEYNLSGTTKKNFSSLDLLNCYQMETVRS
DNIIFQFTKCCPKPKDKSNLLVFRNNGVSDVPTSPQLRPTHMNMVFMHFKIRNEDLIFNESLGGTFTKIF
KGVRRVVDYDGLHETEVLLKVLDAHRNYSSEFFEAASMMSKLSHKHLVNLNYGVCVCGDENILVQEFVK
FGSLDLYLKKKNCINILWKLEVAKQLAWAMHFEENTLIHGNCVCAKNILLIREEDRKTGNPPFIKLSDP
GISITVLPKDILQERIPWVPECIENPKNLNLATDKWSFGTTLWEICSGGDKPLSALDSQRKLQFYEDRH
QLPAPKWAELANLINNCMDYEPDFRPSFRAIRDNLNLSLFTPDYELLTENDMLPNMRIGALGFGS
GAFEDRDP
PTQFEERHLKFLQLGKGNFGSVEMCRYDPLQDNTGEVVAVKKLQHSTEEHLRDFEREIEILKSLQHDNI
VKYKGVCSAGRRNLKIMEYLPYGLRDYLQKHKERIDHIKLLQYTSQICKGMEYLGTKRYIHRDLATR
NILVENENRVKIGDFGLTKVLPQDKEYYKVEKPGESPIFWYAPESLTSKFSVASDVSFGVVLVELFTY
IEKSKSPPAEFMRMIGNDKQGQMIVFHLLIELLKNNGRLLPRPDGCPDEIYMIMTECWNNNVNQRPSFRDL
A
LRVDQIRDNMAG


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	130.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining



[View online »](#)

Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004963
Locus ID:	3717
UniProt ID:	O60674
RefSeq Size:	5097
Cytogenetics:	9p24.1
RefSeq ORF:	3396
Synonyms:	JTK10
Summary:	<p>This gene encodes a non-receptor tyrosine kinase that plays a central role in cytokine and growth factor signalling. The primary isoform of this protein has an N-terminal FERM domain that is required for erythropoietin receptor association, an SH2 domain that binds STAT transcription factors, a pseudokinase domain and a C-terminal tyrosine kinase domain. Cytokine binding induces autophosphorylation and activation of this kinase. This kinase then recruits and phosphorylates signal transducer and activator of transcription (STAT) proteins. Growth factors like TGF-beta 1 also induce phosphorylation and activation of this kinase and translocation of downstream STAT proteins to the nucleus where they influence gene transcription. Mutations in this gene are associated with numerous inflammatory diseases and malignancies. This gene is a downstream target of the pleiotropic cytokine IL6 that is produced by B cells, T cells, dendritic cells and macrophages to produce an immune response or inflammation. Disregulation of the IL6/JAK2/STAT3 signalling pathways produces increased cellular proliferation and myeloproliferative neoplasms of hematopoietic stem cells. A nonsynonymous mutation in the pseudokinase domain of this gene disrupts the domains inhibitory effect and results in constitutive tyrosine phosphorylation activity and hypersensitivity to cytokine signalling. This gene and the IL6/JAK2/STAT3 signalling pathway is a therapeutic target for the treatment of excessive inflammatory responses to viral infections. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Jul 2020]</p>
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Adipocytokine signaling pathway, Chemokine signaling pathway, Jak-STAT signaling pathway

Product images:

Coomassie blue staining of purified JAK2 protein (Cat# TP320503). The protein was produced from HEK293T cells transfected with JAK2 cDNA clone (Cat# [RC220503]) using MegaTran 2.0 (Cat# [TT210002]).