

## Product datasheet for **TP320517M**

### NM23A (NME1) (NM\_198175) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human non-metastatic cells 1, protein (NM23A) expressed in (NME1), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220517 representing NM_198175 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MVLLSTLGIVFQGEPPISSCDTGTMANCERTFIAIKPDGVQRGLVGEIIRFEQKGFRLVGLKFMQASE DLLKEHYVDLKDRPFFAGLVKYMHSGPVVAMVWEGLNVVKTGRVMLGETNPADSKPGTIRGDFCIQVGRN IIHGSDSVESAEKEIGLWFHPEELVDYTSCAQNWIYE  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	19.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_937818</a>
Locus ID:	4830
UniProt ID:	<a href="#">P15531</a>



[View online »](#)

RefSeq Size: 1031

Cytogenetics: 17q21.33

RefSeq ORF: 531


Synonyms: AWD; GAAD; NB; NBS; NDKA; NDPK-A; NDPKA; NM23; NM23-H1

**Summary:** This gene (NME1) was identified because of its reduced mRNA transcript levels in highly metastatic cells. Nucleoside diphosphate kinase (NDK) exists as a hexamer composed of 'A' (encoded by this gene) and 'B' (encoded by NME2) isoforms. Mutations in this gene have been identified in aggressive neuroblastomas. Two transcript variants encoding different isoforms have been found for this gene. Co-transcription of this gene and the neighboring downstream gene (NME2) generates naturally-occurring transcripts (NME1-NME2), which encodes a fusion protein comprised of sequence sharing identity with each individual gene product. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome, Stem cell - Pluripotency

**Protein Pathways:** Metabolic pathways, Purine metabolism, Pyrimidine metabolism

### Product images:



Coomassie blue staining of purified NME1 protein (Cat# [TP320517]). The protein was produced from HEK293T cells transfected with NME1 cDNA clone (Cat# [RC220517]) using MegaTran 2.0 (Cat# [TT210002]).