

Product datasheet for TP320550M

VAV3 (NM_006113) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human vav 3 guanine nucleotide exchange factor (VAV3), transcript variant 1, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220550 representing NM_006113 Red=Cloning site Green=Tags(s)

MEPWKQCAQWLIHCKVLPTNHRVTWDSAQVFDLAQTLRDGVLLCQLLNNLRAHSINLKEINLRPQMSQFL
CLKNIRFTLTACCETFGMRKSELFEAFDLFDVRDFGKVIETLSRLSRTPIALATGIRPFPTESINDEDI
YKGLPDLIDETLVEDEEDLYDCVYGEDEGGEVYEDLMKAEAAHQPKCPENDIRSCCLAEIKQTEEKYTET
LESIEKYFMAPLKRFLTAAEFDSVFINIPELVKLRNLMQEIHDSIVNKNDQNLQVFINYKERLVIYQG
YCSGVESAISLDYISKTKEDVKLLEEC SKRANNGKFTLRDLLVPMQRVLKYHLLQLV KHHTDPTE
KANLKLALDAMKDLAQYVNEVKRDNETLREIKQFQLSIENLNQPVLLFGRPQGDGEIRITTLDKHTKQER
HIFLFDLAVIVCKRKGDNEMKEIIDLQQYKIANNPTTDKENKKWSYGFYLIHTQGQNGLEFYCKTKDLK
KKWLEQFEMALS NIRPDYADSNFHFDFKMHTFTRVTSCKVCQMLLRGTFYQGYLCFKCGARAHKECLGRVD
NCGRVNSGEQGTLLKPEKRTNGLRRTPKQVDPGLPKMQVIRNYSGTPPALHEGPSLHLQAGDTELLKG
DAHSLFWQGRNLASGEVGFPSDAVKPCPCVPKVDYSCQPWYAGAMERLQAETELINRVNSTYLVRHRT
KESGEYAIKYNNEAKHIKILTRDGFHIAENRKFSLMELVEYKHHSLKEGFRITLDTTLQFPYKEPE
HSAGQRGNRAGNSLLSPKVLGIAIARYDFCARDMRELSLLKGDVVKIYTKMSANGWWRGEVNGRVGWFPS
TYVEEDE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	97.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006104
Locus ID:	10451
UniProt ID:	Q9UKW4
RefSeq Size:	4768
Cytogenetics:	1p13.3
RefSeq ORF:	2541
Summary:	This gene is a member of the VAV gene family. The VAV proteins are guanine nucleotide exchange factors (GEFs) for Rho family GTPases that activate pathways leading to actin cytoskeletal rearrangements and transcriptional alterations. This gene product acts as a GEF preferentially for RhoG, RhoA, and to a lesser extent, RAC1, and it associates maximally with the nucleotide-free states of these GTPases. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome
Protein Pathways:	B cell receptor signaling pathway, Chemokine signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Leukocyte transendothelial migration, Natural killer cell mediated cytotoxicity, Regulation of actin cytoskeleton, T cell receptor signaling pathway

Product images:



Coomassie blue staining of purified VAV3 protein (Cat# [TP320550]). The protein was produced from HEK293T cells transfected with VAV3 cDNA clone (Cat# [RC220550]) using MegaTran 2.0 (Cat# [TT210002]).