

Product datasheet for TP320566L

CTPS2 (NM_175859) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human CTP synthase II (CTPS2), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220566 protein sequence Red=Cloning site Green=Tags(s)

MKYILVTGGVISGIGKGIASSIGTILKSCGLRVTAIKIDPYINIDAGTFSPYEHGEVFLNDGGGEVDLD
LGNYERFLDINLYKDNNITGKIYQHVINKERRGDYLGKTVQVPHITDAVQEWMNQAKVPVDGNKEEP
QICVIELGGTIGDIEGMPFVEAFRQFQFKAKRENFNCNIHVSLVPLSATGEQKTKPTQNSVRALRGLGLS
PDLIVCRSSTPIEMAVKEKISMFCNVNPEQVICIHDVSSTYRVPVLEEVSIVKYFKERLHLPIDGSASN
LLFKWRNMADRYERLQKICSIALVGGYTKLRDCYASVFKALEHSALAINHKLNLMYIDSIDLEKITETED
PVKFHEAWQKLCCKADGILVPGGFGIRGTLGKLQAISWARTKKIPFLGVCLGMQLAVIEFARNCLNLKDAD
STEFRPNAPVPLVIDMPEHNPGNLGGTMRLGIRRTVFKTENSILRKLYGDVPPFIEERHRHRFEVNPNLK
QFEQNDLSFVGQDQDVGDRMEIILANHPYFVGQVQFHPEFSSRPMKPSPPYLGLLLAATGNLNAYLQQGCK
LSSSDRYSDASDDSFSEPRIAELEIS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	65.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_787055
Locus ID:	56474
UniProt ID:	Q9NRF8 , A0A024RC00
RefSeq Size:	3973
Cytogenetics:	Xp22.2
RefSeq ORF:	1758
Synonyms:	GATD5B
Summary:	The protein encoded by this gene catalyzes the formation of CTP from UTP with the concomitant deamination of glutamine to glutamate. This protein is the rate-limiting enzyme in the synthesis of cytosine nucleotides, which play an important role in various metabolic processes and provide the precursors necessary for the synthesis of RNA and DNA. Cancer cells that exhibit increased cell proliferation also exhibit an increased activity of this encoded protein. Thus, this protein is an attractive target for selective chemotherapy. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]
Protein Pathways:	Metabolic pathways, Pyrimidine metabolism

Product images:



Coomassie blue staining of purified CTPS2 protein (Cat# [TP320566]). The protein was produced from HEK293T cells transfected with CTPS2 cDNA clone (Cat# [RC220566]) using MegaTran 2.0 (Cat# [TT210002]).