

## Product datasheet for **TP320589**

### ALR (GFER) (NM\_005262) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human growth factor, augments liver regeneration (GFER), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC220589 representing NM_005262
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MAAPGERGRFHGGNLFPLPGGARSEMDDLATDARGRGAGRRDAAASASTPAQAPTS DSPVAEDASRRRP  
CRACVDFKTWMRTQQRDTKRFREDCPPDREELGRHSWAVLHTLAAYPDLPTPEQQQDMAQFIHLFSKFY  
PCEECAEDLRKRLCRNHPDTRTRACFTQWLCHLHNEVNRKLGKPDFDCSKVDERWRD GWKDGSCD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	23.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_005253</a>
Locus ID:	2671
UniProt ID:	<a href="#">P55789</a>
RefSeq Size:	2447



[View online »](#)


Cytogenetics: 16p13.3

RefSeq ORF: 615

Synonyms: ALR; ERV1; HERV1; HPO; HPO1; HPO2; HSS; MMCHD; MPMCD

**Summary:** The hepatotrophic factor designated augments liver regeneration (ALR) is thought to be one of the factors responsible for the extraordinary regenerative capacity of mammalian liver. It has also been called hepatic regenerative stimulation substance (HSS). The gene resides on chromosome 16 in the interval containing the locus for polycystic kidney disease (PKD1). The putative gene product is 42% similar to the scERV1 protein of yeast. The yeast scERV1 gene had been found to be essential for oxidative phosphorylation, the maintenance of mitochondrial genomes, and the cell division cycle. The human gene is both the structural and functional homolog of the yeast scERV1 gene. [provided by RefSeq, Jul 2008]

### Product images:



Coomassie blue staining of purified GFER protein (Cat# TP320589). The protein was produced from HEK293T cells transfected with GFER cDNA clone (Cat# [RC220589]) using MegaTran 2.0 (Cat# [TT210002]).