

## Product datasheet for TP320599

### PKC beta 1 (PRKCB) (NM\_212535) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protein kinase C, beta (PRKCB), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220599 representing NM_212535 Red=Cloning site Green=Tags(s)

MADPAAGPPPSEGEESTVRFARKGALRQKNVHEVKNHKFTARFFKQPTFCSHCTDFIWGFGKQGFQCQV  
C  
CFVWHKRCHEFVTFSCPGADKGPASDDPRSKHKFKIHTYSSPTFCDHCGSLLYGLIHQGMKCDTCMMNV  
H  
KRCVMNVPSLCGTDHTERRGRIYIQAHDRLVIVLVRDAKNLVPMDPNGLSDPYVKLKLIPDPKSESKQ  
KTKTIKCSLNPEWNETFRFQLKESDKDRRLSVEIWDWDLTSRNDFMGSLSFGISELQKASVDGWFKLLSQ  
EEGEYFNVPVPEGSEANEELRQKFERAKISQGTKVPEEKTTNTVSKFDNNGNRDRMKLTDNFNLMVLGK  
GSFGKVMLSERKGTDELYAVKILKDVVIQDDDVECTMVEKRVLALPGKPPFLTQLHSCFQTMDRLYFVM  
EYVNGGDLMYHIQQVGRFKEPHAVFYAAEIAIGLFFLQSKGIYRDLKLDNVMLDSEGHIKIADFGMCKE  
NIWDGVTTKTCGTPDYIAPEIIAYQPYGKSVDWWAFGVLLYEMLAGQAPFEGEDEDELQFSIMEHNVAY  
PKSMSKEAVAICKGLMTKHPGKRLGCGPEGERDIKEHAFFRYIDWEKLERKEIQPPYKPKARDKRDSNF  
DKEFTRQPVELTPTDKLFIMNLDQNEFAGFSYTNPEFVINV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	76.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_997700</a>
<b>Locus ID:</b>	5579
<b>UniProt ID:</b>	<a href="#">P05771</a>
<b>RefSeq Size:</b>	2319
<b>Cytogenetics:</b>	16p12.2-p12.1
<b>RefSeq ORF:</b>	2013
<b>Synonyms:</b>	PKC-beta; PKCB; PKCbeta; PKCI(2); PRKCB1; PRKCB2
<b>Summary:</b>	<p>Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways. PKC family members also serve as major receptors for phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play a distinct role in cells. The protein encoded by this gene is one of the PKC family members. This protein kinase has been reported to be involved in many different cellular functions, such as B cell activation, apoptosis induction, endothelial cell proliferation, and intestinal sugar absorption. Studies in mice also suggest that this kinase may also regulate neuronal functions and correlate fear-induced conflict behavior after stress. Alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Jul 2008]</p>
<b>Protein Families:</b>	Druggable Genome, Protein Kinase
<b>Protein Pathways:</b>	B cell receptor signaling pathway, Calcium signaling pathway, Chemokine signaling pathway, ErbB signaling pathway, Fc epsilon RI signaling pathway, Fc gamma R-mediated phagocytosis, Focal adhesion, Gap junction, Glioma, GnRH signaling pathway, Leukocyte transendothelial migration, Long-term depression, Long-term potentiation, MAPK signaling pathway, Melanogenesis, Natural killer cell mediated cytotoxicity, Non-small cell lung cancer, Pathways in cancer, Phosphatidylinositol signaling system, Tight junction, Vascular smooth muscle contraction, VEGF signaling pathway, Vibrio cholerae infection, Wnt signaling pathway

**Product images:**

Coomassie blue staining of purified PRKCB protein (Cat# TP320599). The protein was produced from HEK293T cells transfected with PRKCB cDNA clone (Cat# [RC220599]) using MegaTran 2.0 (Cat# [TT210002]).