

Product datasheet for **TP320657L**

C10orf53 (NM_001042427) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 10 open reading frame 53 (C10orf53), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220657 representing NM_001042427 Red =Cloning site Green =Tags(s) MPKNAVILRYGPYSAAGLPVEHHTFRLQGLQAVLAIDGHEVILEKIEDWNVWELMVNEEVIFHCNIKDL EFGGDGKLDPLCEKARIAVLNAY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	10.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001035892
Locus ID:	282966
UniProt ID:	Q8N6V4
RefSeq Size:	642




[View online »](#)

Cytogenetics: 10q11.23

RefSeq ORF: 279

Product images:



Coomassie blue staining of purified C10orf53 protein (Cat# [TP320657]). The protein was produced from HEK293T cells transfected with C10orf53 cDNA clone (Cat# [RC220657]) using MegaTran 2.0 (Cat# [TT210002]).