

Product datasheet for TP320667

RIOK1 (NM_031480) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RIO kinase 1 (yeast) (RIOK1), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220667 representing NM_031480 Red =Cloning site Green =Tags(s)

MDYRLLMSRVPGQFDDADSSDSENRLKTVKEKDDILFEDLQDNVNENGEIEDEEEEGYDDDDDDW
DWDEGVGKLAKGYVWNGGSPQANRQTSDDSSAKMSTPADKVLRFKFNKINLDKLNVTDSVINKVTEKSR
QKEADMYRIKDKADRATVEQVLDPRTRMILFKMLTRGIITEINGCISTGKEANVYHASTANGESRAIKY
KTSILVFKDRDKYVSGEFRFRHGYCKGNPRKMVKTWAEKEMRNLRLNTAEIPCPEPIMLRSHVLVMSFI
GKDDMPAPLLKNVQLSESKARELYLQVIQYMRRMYQDARLVHADLSEFNMLYHGGGVYIIDVSQSVEHDH
PHALEFLRKDCANVNDFFMRHSVAVMTVRELFVFTDPSITHENMDAYLSKAMEIASQRTKEERSSQDHV
DEEVFKRAYIPRTLNEVKNYERDMDIIMKLKEEDMAMNAQQDNILYQTVTGLKKDLSGVQKVPALLENQV
EERTCSDSEDIGSSECSDDTSEEQGDHARPKKHTTDPDIDKKERKMMVKEAQREKRKNKIPKHVKKRKEK
TAKTKKGGK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	65.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_113668
Locus ID:	83732
UniProt ID:	Q9BRS2
RefSeq Size:	2495
Cytogenetics:	6p24.3
RefSeq ORF:	1704
Synonyms:	AD034; bA288G3.1; RIO1; RRP10
Summary:	The protein encoded by this gene competes with pICln for inclusion in the protein arginine methyltransferase 5 complex. This complex targets substrates for dimethylation. The encoded protein is essential for the last steps in the maturation of 40S subunits. [provided by RefSeq, Jan 2017]
Protein Families:	Druggable Genome, Protein Kinase

Product images:



Coomassie blue staining of purified RIOK1 protein (Cat# TP320667). The protein was produced from HEK293T cells transfected with RIOK1 cDNA clone (Cat# [RC220667]) using MegaTran 2.0 (Cat# [TT210002]).