

Product datasheet for TP320718

C20orf112 (NOL4L) (NM_080616) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens chromosome 20 open reading frame 112 (C20orf112), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC220718 representing NM_080616
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MSDSTWMSADPHLASSLSPSQDERMRSPQNLHSQEDDDSSSESGSGNGSSTLNPSTSSSTQGDPAFPEMN
GNGAVAPMDFTTAAEDQPINLCKLPPATALGTASYPSDGCGADGLRSRVKYGVKTTPEPPYSSGSYDS
IKTEVSGPEDLTVGRAPTADDDDDHDDHEDNDKMNDSEGMDPERLKAFNMFVRLFVDENLDRMVPISK
QPKEEIQAIIESCQRQFPEFQERARKRIRTYLKSCRRMKKNGMEMTRPTPPHLTSAMAENILAAACESET
RKAARKMRLEIYQSSQDEPIALDKQHSRDSAAITHSTYSLPASSYSQDPVYANGGLNYSYRGYGALSSNL
QPPASLQGTGNHNGPTDLSMKGGASTTSTTPTPTPSSTSTSRPVPTAQLSPTEISAVRQLIAGYRESAAF
LLRSADELENLILQQN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	47 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_542183](#)

Locus ID: 140688

UniProt ID: [Q96MY1](#)


RefSeq Size: 2100

Cytogenetics: 20q11.21

RefSeq ORF: 1308

Synonyms: C20orf112; C20orf113

Product images:



Coomassie blue staining of purified NOL4L protein (Cat# TP320718). The protein was produced from HEK293T cells transfected with NOL4L cDNA clone (Cat# [RC220718]) using MegaTran 2.0 (Cat# [TT210002]).