

Product datasheet for **TP320737L**

KIAA0664 (CLUH) (NM_015229) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens KIAA0664 (KIAA0664), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220737 representing NM_015229 Red =Cloning site Green =Tags(s)

MLLNGDCPESLKKEAAAAEPPRENLDEAGPGDETTGQEVIVIQDTGFSVKILAPGIEPFSLQVSPQEMV
 QEIHQVLMRDREDTCHRTCSLHLDGNVLDHFSELRSVEGLQEGSVLRVVEEPTVREARIHVHRVDRLLK
 SLDPDAFNGVDCNSLSFLSVFTDGDLDGDSGKRKKGLEMDPIDCTPPEYILPGSRERPLCPLQPQNRDWK
 PLQCLKVLTMSGWNPPPGNRKMHGDLMYLFVITAEDRQVSITASTRGFYLNQSTAYHFNPKPASPRFLSH
 SLVELLNQISPTFKKNFAVLQKKRVQRHPFERIATPFQVYSWTAPQAEHAMDCVRAEDAYTSRLGYEEHI
 PGQTRDWNHEELQTTRELPRKNLPERLLRERAIKVFHSDFTAAATRGAMAVIDGNVMAINPSEETKMQMFI
 WNNIFFSLGFDVRDHYKDFGGDVAAYVAPTNDLNGVRTYNAVDVEGLYLTGTVVVDYRGYRVTAQSIIPG
 ILERDQEQSVIYGSIDFGKTVVSHPRYLELLERTSRPLKILRHQVLNDRDEEVELCSSVECKGIIGNDGR
 HYILDLLRTFPPDLNFLPVPGEELPEECARAGFPRAHRHKLCCLRQELVDAFVEHRYLLFMKLAALQLMQ
 QNASQLETPSSLENGGPSSLESKSEDPGQGEAGSEEEGSSASGLAKVKELAETIAADDGTDPRSREVIRN
 ACKAVGSISSTAFDIRFNPDIFSPGVRFPESCQDEV RDQKQLLKDA AFLSCQIPGLVKDCMEHAVLPV
 DGATLAEVMRQRGINMRYLGKVL ELVLRSPARHQLDHVF KIGIGELITRSAKHIFKTYLQGVLSGLSAA
 ISHFLNCF LSSYPNPVAHLPADELVSKKRNRKRNRPPGAADNTAWAVMTPQELWKNICQEAKNYFD F D
 L
 ECETVDQAVETYGLQKITLLREISLKTGIQVLLKEYSFD SRHKPAFTEEDVLNIFPVVKHVNPKASDAFH
 FFQSGQAKVQQGFLKEGCELINEALNLFNNVYGAMHVETCACLRLRLHYIMGDYAEALS NQQKAVLM
 S
 ERVMGTEHPNTIQEYMH LALYCFASSQLSTALSLLYRARYLM LLVFGEDHPEMALLDNNIGLV LHGVMEY
 DLSLRFLENALAVSTKYHGP KALKVALSHHLVARVYESKA EFRSALQHEKEGYTIYKTQLGEDHEKTKE S
 SEYLKCLTQQAVALQRTMNEIYRNGSSANIPPLKFTAPSMASVLEQLNVINGILFIPLSQKDLENLKA EV
 ARRHQLQEASRNDRAE EPMATEPAPAGAPGDLGSQPPAAKDPSPSVQGG

SGPTRRRLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	146.5 kDa



[View online »](#)

Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_056044
Locus ID:	23277
UniProt ID:	O75153
RefSeq Size:	5252
Cytogenetics:	17p13.3
RefSeq ORF:	3927
Synonyms:	CLU1
Summary:	mRNA-binding protein involved in proper cytoplasmic distribution of mitochondria. Specifically binds mRNAs of nuclear-encoded mitochondrial proteins in the cytoplasm and regulates transport or translation of these transcripts close to mitochondria, playing a role in mitochondrial biogenesis.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified CLUH protein (Cat# [TP320737]). The protein was produced from HEK293T cells transfected with CLUH cDNA clone (Cat# [RC220737]) using MegaTran 2.0 (Cat# [TT210002]).