

## Product datasheet for **TP320853**

### DHRS9 (NM\_199204) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human dehydrogenase/reductase (SDR family) member 9 (DHRS9), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220853 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MLFWVLGLLILCGFLWTRKGLKIEDITDKYIFITGCDSGFGNLAARTFDKKGPHVIAACLTESGSTALK  
AETSERLRTVLLDVTDPENVKRTAQWVKVQVGEKGLWGLINNAGVPGVLAPTDWLTLEDYREPIEVNLFGL  
LISVTLNMLPLVKKKAQGRVINVSSVGGRLAIVGGGYTPSKYAVEGFNDSLRRDMKAFGVHVSCIEPGLFK  
TNLADPVKVIKLAWEQLSPDIKQQYGEYIEKSLDKLKGKNSYVNMDLSPVVECMDHALTSLFPKTH  
YAAGKDAKIFWIPLSHMPAALQDFLLKQKAELANPKAV

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**


Tag:	C-Myc/DDK
Predicted MW:	35 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_954674</a>
Locus ID:	10170



[View online »](#)

UniProt ID:	<a href="#">Q9BPW9</a>
RefSeq Size:	1993
Cytogenetics:	2q31.1
RefSeq ORF:	957
Synonyms:	3-alpha-HSD; 3ALPHA-HSD; RDH-E2; RDH-TBE; RDH15; RDHL; RDHTBE; RETSDR8; SDR9C4
Summary:	This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family. The encoded protein has been identified as a moonlighting protein based on its ability to perform mechanistically distinct functions. This protein demonstrates oxidoreductase activity toward hydroxysteroids and is able to convert 3-alpha-tetrahydroprogesterone to dihydroxyprogesterone and 3-alpha-androstanediol to dihydroxyprogesterone in the cytoplasm, and may additionally function as a transcriptional repressor in the nucleus. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Retinol metabolism

### Product images:



Coomassie blue staining of purified DHRS9 protein (Cat# TP320853). The protein was produced from HEK293T cells transfected with DHRS9 cDNA clone (Cat# [RC220853]) using MegaTran 2.0 (Cat# [TT210002]).