

Product datasheet for TP320981

Nogo A (RTN4) (NM_207520) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human reticulon 4 (RTN4), transcript variant 4, 20 µg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC220981 representing NM_207520
Red=Cloning site Green=Tags(s)

MEDLDQSPLVSSSDSPRPQPAFKYQFVREPEDEEEEEEEEEDEDEDLEELEVLERKPAAGLSAAPVPT
 APAAGAPLMDFGNDVFPAPRGPLPAAPPVAPERQPSWDPSVSTVPAPSPLSAAAVSPSKLPEDDEPP
 ARPPPPPPASVSPAEPVWTPPAPAPAAPPSTPAAPKRRGSSGSVDETLFALPAASEPVIRSSAVDLLY
 WRDIKKTGVVFGASLFLLSLTVFSIVSVTAYIALALLSVTISFRIYKGMVQAIQKSDEGHPFRAYLESE
 VAISEELVQKYSNSALGHVNCTIKELRRLFLVDDLVDLSLKFAVLMWVFTYVGFNGLTLILALISLFS
 VPVIYERHQAQIDHYLGLANKNVKDAMAKIQAKIPGLKRKAE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 42.1 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage: Store at -80°C.
Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq: [NP_997403](#)
Locus ID: 57142




[View online »](#)

UniProt ID:	Q9NQC3
RefSeq Size:	2471
Cytogenetics:	2p16.1
RefSeq ORF:	1176
Synonyms:	ASY; Nbla00271; Nbla10545; NI220/250; NOGO; NSP; NSP-CL; RTN-X; RTN4-A; RTN4-B1; RTN4-B2; RTN4-C

Summary: This gene belongs to the family of reticulon encoding genes. Reticulons are associated with the endoplasmic reticulum, and are involved in neuroendocrine secretion or in membrane trafficking in neuroendocrine cells. The product of this gene is a potent neurite outgrowth inhibitor which may also help block the regeneration of the central nervous system in higher vertebrates. Alternatively spliced transcript variants derived both from differential splicing and differential promoter usage and encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified RTN4 protein (Cat# TP320981). The protein was produced from HEK293T cells transfected with RTN4 cDNA clone (Cat# [RC220981]) using MegaTran 2.0 (Cat# [TT210002]).