

Product datasheet for TP320989M

PKN2 (NM_006256) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protein kinase N2 (PKN2), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC220989 representing NM_006256 Red=Cloning site Green=Tags(s)

MASNPERGEILLTELQGDSRSLPFSENVSAVQKLDIFSMTMVQQKLDLDIKDRIKREIRKELKIKEGAENLR
KVTTDKKSLAYVDNLIKKSNNKLEELHHLQELNAHIVSDPEDITDCPRTPDTPNNDPRCSTSNRLKA
LQKQLDIELKVKQGAENMIQMYSNGSSKDRKLHGTAQQLLQDSKTKIEVIRMQILQAVQTNELAFDPAK
VISPLELRMEELRHHFRIEFAVAEGAKNVMKLLGSGKVTDRKALSEAQARFNESSQKDLLKYSLEQRLN
EVPKNHPKSRIIIIEELSLVAASPTLSPRQSMISTQNQYSTLSKPAALTGTLEVRMLMGCQDILENVPGRSK
ATSVALPGWSPSETRSSFMSRTSKSKSGSSRNLKTDDLSDNVCVAVLKLNTVVGQTSWKPISNQSWDQK
FTLELDRSRELEISVYWRDWRSLCAVKFLRLEDFLDNQRHGMCLYLEPQGTLFAEVTFFNPVIERRPKLQ
RQKKIFSKQQGKTLRAPQMNINIATWGRVLRRAIPTVNHSGTFSPQAPVPTTVPVVDVRIPLAPPASD
STVTKLDFDLEPEPPPAPPRASSLGEIDESSELRVLDIPGQDSETVFDIQNDRNSILPKSQSEYKPDTPQ
SGLEYSIGIQELEDRRSQRFQFNLQDFRCCAVLGRGHFGKVLAEYKNTNEMFAIKALKKGDIVARDEVD
SLMCEKRIFETVNSVRHPFLVNLFACFQTKHEVCFVMEYAAGGDLMMHIHTDVFSEPRAVFYAACVVLGL
QYLHEHKIVYRDCLKLDNLLDTEGFVKIADFLGCKEGMGYGDRTSTFCGTPEFLAPEVLTETSYTRAVDW
WGLGVLIYEMLVGESPFPGDDEEEVFDSIVNDEVRYRFLSTEAISIMRLLRRNPERRLGASEKDAEDV
KKHPFFRLIDWSALMDKKVPPFIPTIRGREDVSNFDDEFTSEAPILTPPREPRILSEEEQEMFRDFDYI
ADWC


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	111.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006247
Locus ID:	5586
UniProt ID:	Q16513
RefSeq Size:	3255
Cytogenetics:	1p22.2
RefSeq ORF:	2952
Synonyms:	Pak-2; PAK2; PRK2; PRKCL2; PRO2042; STK7
Summary:	<p>PKC-related serine/threonine-protein kinase and Rho/Rac effector protein that participates in specific signal transduction responses in the cell. Plays a role in the regulation of cell cycle progression, actin cytoskeleton assembly, cell migration, cell adhesion, tumor cell invasion and transcription activation signaling processes. Phosphorylates CTTN in hyaluronan-induced astrocytes and hence decreases CTTN ability to associate with filamentous actin. Phosphorylates HDAC5, therefore lead to impair HDAC5 import. Direct RhoA target required for the regulation of the maturation of primordial junctions into apical junction formation in bronchial epithelial cells. Required for G2/M phases of the cell cycle progression and abscission during cytokinesis in a ECT2-dependent manner. Stimulates FYN kinase activity that is required for establishment of skin cell-cell adhesion during keratinocytes differentiation. Regulates epithelial bladder cells speed and direction of movement during cell migration and tumor cell invasion. Inhibits Akt pro-survival-induced kinase activity. Mediates Rho protein-induced transcriptional activation via the c-fos serum response factor (SRF). Involved in the negative regulation of ciliogenesis (PubMed:27104747).[UniProtKB/Swiss-Prot Function]</p>
Protein Families:	Druggable Genome, Protein Kinase

Product images:

Coomassie blue staining of purified PKN2 protein (Cat# [TP320989]). The protein was produced from HEK293T cells transfected with PKN2 cDNA clone (Cat# [RC220989]) using MegaTran 2.0 (Cat# [TT210002]).