

# **Product datasheet for TP321018**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## NUMB (NM\_001005744) Human Recombinant Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human numb homolog (Drosophila) (NUMB), transcript variant 2, 20

μg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC221018 representing NM\_001005744

or AA Sequence: Red=Cloning site Green=Tags(s)

MNKLRQSFRRKKDVYVPEASRPHQWQTDEEGVRTGKCSFPVKYLGHVEVDESRGMHICEDAVKRLKAERK FFKGFFGKTGKKAVKAVLWVSADGLRVVDEKTKDLIVDQTIEKVSFCAPDRNFDRAFSYICRDGTTRRWI CHCFMAVKDTGERLSHAVGCAFAACLERKQKREKECGVTATFDASRTTFTREGSFRVTTATEQAEREEIM KQMQDAKKAETDKIVVGSSVAPGNTAPSPSSPTSPTSDATTSLEMNNPHAIPRRHAPIEQLARQGSFRGF PALSQKMSPFKRQLSLRINELPSTMQRKTDFPIKNAVPEVEGEAESISSLCSQITNAFSTPEDPFSSAPM TKPVTVVAPQSPTFQGTEWGQSSGAASPGLFQAGHRRTPSEADRWLEEVSKSVRAQQPQASAAPLQPVLQ PPPPTAISQPASPFQGNAFLTSQPVPVGVVPALQPAFVPAQSYPVANGMPYPAPNVPVVGITPSQMVANV FGTAGHPQAAHPHQSPSLVRQQTFPHYEASSATTSPFFKPPAQHLNGSAAFNGVDDGRLASADRHTEVPT

GTCPVDPFEAQWAALENKSKQRTNPSPTNPFSSDLQKTFEIEL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

**Predicted MW:** 65.7 kDa

**Concentration:** >0.05 μg/μL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.





#### NUMB (NM\_001005744) Human Recombinant Protein - TP321018

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

**RefSeq:** NP 001005744

**Locus ID:** 8650

UniProt ID: <u>P49757</u>, <u>A0A024R681</u>

RefSeq Size: 3503

**Cytogenetics:** 14q24.2-q24.3

RefSeq ORF: 1809

**Synonyms:** C14orf41; c14\_5527; S171

**Summary:** The protein encoded by this gene plays a role in the determination of cell fates during


development. The encoded protein, whose degradation is induced in a proteasome-dependent

manner by MDM2, is a membrane-bound protein that has been shown to associate with EPS15, LNX1, and NOTCH1. Alternative splicing results in multiple transcript variants. [provided

by RefSeq, Feb 2016]

**Protein Pathways:** Notch signaling pathway

## **Product images:**



Coomassie blue staining of purified NUMB protein (Cat# TP321018). The protein was produced from HEK293T cells transfected with NUMB cDNA clone (Cat# [RC221018]) using

MegaTran 2.0 (Cat# [TT210002]).