

## **Product datasheet for TP321022**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

#### DKK3 (NM 015881) Human Recombinant Protein

**Product data:** 

**Product Type:** Recombinant Proteins

**Description:** Recombinant protein of human dickkopf homolog 3 (Xenopus laevis) (DKK3), transcript variant

1, 20 μg

Species: Human Expression Host: HEK293T

**Expression cDNA** >RC221022 protein sequence **Clone or AA Sequence:** Red=Cloning site Green=Tags(s)

MQRLGATLLCLLLAAAVPTAPAPAPTATSAPVKPGPALSYPQEEATLNEMFREVEELMEDTQHKLRSAVE EMEAEEAAAKASSEVNLANLPPSYHNETNTDTKVGNNTIHVHREIHKITNNQTGQMVFSETVITSVGDEE GRRSHECIIDEDCGPSMYCQFASFQYTCQPCRGQRMLCTRDSECCGDQLCVWGHCTKMATRGSNGTICDN QRDCQPGLCCAFQRGLLFPVCTPLPVEGELCHDPASRLLDLITWELEPDGALDRCPCASGLLCQPHSHSL VYVCKPTFVGSRDQDGEILLPREVPDEYEVGSFMEEVRQELEDLERSLTEEMALGEPAAAAAALLGGEEI

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV** 

Tag: C-Myc/DDK

**Predicted MW:** 36.2 kDa

**Concentration:** >0.05 µg/µL as determined by microplate BCA method

**Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining

**Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

**Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

**Note:** For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling

conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 056965

Locus ID: 27122



## DKK3 (NM\_015881) Human Recombinant Protein - TP321022

UniProt ID: Q9UBP4

RefSeq Size: 2769

Cytogenetics: 11p15.3 RefSeq ORF: 1050


Synonyms: REIC; RIG

**Summary:** This gene encodes a protein that is a member of the dickkopf family. The secreted protein

contains two cysteine rich regions and is involved in embryonic development through its interactions with the Wnt signaling pathway. The expression of this gene is decreased in a variety of cancer cell lines and it may function as a tumor suppressor gene. Alternative splicing results in multiple transcript variants encoding the same protein. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome, Secreted Protein

# **Product images:**



Coomassie blue staining of purified DKK3 protein (Cat# TP321022). The protein was produced from HEK293T cells transfected with DKK3 cDNA clone (Cat# [RC221022]) using MegaTran 2.0 (Cat# [TT210002]).