


Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_683706
Locus ID:	114883
UniProt ID:	Q96SU4 , Q8TAS8
RefSeq Size:	2880
Cytogenetics:	1p32.3
RefSeq ORF:	2169
Synonyms:	ORP-9; ORP9
Summary:	This gene encodes a member of the oxysterol-binding protein (OSBP) family, a group of intracellular lipid receptors. Most members contain an N-terminal pleckstrin homology domain and a highly conserved C-terminal OSBP-like sterol-binding domain, although some members contain only the sterol-binding domain. This family member functions as a cholesterol transfer protein that regulates Golgi structure and function. Multiple transcript variants, most of which encode distinct isoforms, have been identified. Related pseudogenes have been identified on chromosomes 3, 11 and 12. [provided by RefSeq, Jul 2010]

Product images:



Coomassie blue staining of purified OSBPL9 protein (Cat# TP321053). The protein was produced from HEK293T cells transfected with OSBPL9 cDNA clone (Cat# [RC221053]) using MegaTran 2.0 (Cat# [TT210002]).