

Product datasheet for TP321104

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CNN2 (NM_004368) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human calponin 2 (CNN2), transcript variant 1, 20 μg

Species: Human
Expression Host: HEK293T

Expression cDNA >RC221104 representing NM_004368
Clone or AA Red=Cloning site Green=Tags(s)

Sequence:

MSSTQFNKGPSYGLSAEVKNRLLSKYDPQKEAELRTWIEGLTGLSIGPDFQKGLKDGTILCTLMNKLQPG SVPKINRSMQNWHQLENLSNFIKAMVSYGMNPVDLFEANDLFESGNMTQVQVSLLALAGKAKTKGLQSGV DIGVKYSEKQERNFDDATMKAGQCVIGLQMGTNKCASQSGMTAYGTRRHLYDPKNHILPPMDHSTISLQM GTNKCASQVGMTAPGTRRHIYDTKLGTDKCDNSSMSLQMGYTQGANQSGQVFGLGRQIYDPKYCPQGTVA

DGAPSGTGDCPDPGEVPEYPPYYQEEAGY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 33.5 kDa

Concentration: >0.05 μg/μL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling

conditions. Avoid repeated freeze-thaw cycles.

RefSeq: <u>NP 004359</u>

Locus ID: 1265





CNN2 (NM_004368) Human Recombinant Protein - TP321104


UniProt ID: Q99439
RefSeq Size: 2478
Cytogenetics: 19p13.3
RefSeq ORF: 927

Summary: The protein encoded by this gene, which can bind actin, calmodulin, troponin C, and

tropomyosin, may function in the structural organization of actin filaments. The encoded protein could play a role in smooth muscle contraction and cell adhesion. Several pseudogenes of this gene have been identified, and are present on chromosomes 1, 2, 3, 6, 9, 11, 13, 15, 16, 21 and 22. Alternative splicing results in multiple transcript variants encoding different isoforms.

[provided by RefSeq, Jan 2015]

Product images:



Coomassie blue staining of purified CNN2 protein (Cat# TP321104). The protein was produced from HEK293T cells transfected with CNN2 cDNA clone (Cat# [RC221104]) using MegaTran 2.0 (Cat# [TT210002]).