

# Product datasheet for TP321225

## LIPF (NM\_004190) Human Recombinant Protein

### **Product data:**

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn


Product data:	
Product Type:	Recombinant Proteins
Description:	Recombinant protein of human lipase, gastric (LIPF), 20 μg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221225 representing NM_004190 Red=Cloning site Green=Tags(s)
	MWLLLTMASLISVLGTTHGLFGKLHPGSPEVTMNISQMITYWGYPNEEYEVVTEDGYILEVNRIPYGKKN SGNTGQRPVVFLQHGLLASATNWISNLPNNSLAFILADAGYDVWLGNSRGNTWARRNLYYSPDSVEFWAF SFDEMAKYDLPATIDFIVKKTGQKQLHYVGHSQGTTIGFIAFSTNPSLAKRIKTFYALAPVATVKYTKSL INKLRFVPQSLFKFIFGDKIFYPHNFFDQFLATEVCSREMLNLLCSNALFIICGFDSKNFNTSRLDVYLS HNPAGTSVQNMFHWTQAVKSGKFQAYDWGSPVQNRMHYDQSQPPYYNVTAMNVPIAVWNGGKDLLADPQD VGLLLPKLPNLIYHKEIPFYNHLDFIWAMDAPQEVYNDIVSMISEDKK
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	43.2 kDa
Concentration:	>0.05 μg/μL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP 004181</u>
Locus ID:	8513



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

	LIPF (NM_004190) Human Recombinant Protein – TP321225
UniProt ID:	<u>P07098</u>
RefSeq Size:	1365
Cytogenetics:	10q23.31
RefSeq ORF:	1194
Synonyms:	GL; HGL; HLAL
Summary:	This gene encodes gastric lipase, an enzyme involved in the digestion of dietary triglycerides in the gastrointestinal tract, and responsible for 30% of fat digestion processes occurring in human. It is secreted by gastric chief cells in the fundic mucosa of the stomach, and it hydrolyzes the ester bonds of triglycerides under acidic pH conditions. The gene is a member of a conserved gene family of lipases that play distinct roles in neutral lipid metabolism. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2010]
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways	Glycerolipid metabolism, Metabolic pathways

#### **Product images:**



Coomassie blue staining of purified LIPF protein (Cat# TP321225). The protein was produced from HEK293T cells transfected with LIPF cDNA clone (Cat# [RC221225]) using MegaTran 2.0 (Cat# [TT210002]).

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2023 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US