

Product datasheet for TP321328

NOVA1 (NM_006489) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human neuro-oncological ventral antigen 1 (NOVA1), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221328 representing NM_006489 Red =Cloning site Green =Tags(s)
	<p>MMAAAPIQQNGTHTGVPIDLDPPDSRKRPLEAPPEAGSTKRTNTGEDGQYFLKVLIPSYAAGSIIGKGGQ TIVQLQKETGATIKLSKSKDFYPGTTTERVCLIQGTVEALNAVHGFIAEKIREMPQNVAKTEPVSILQPQT TVNPDRIKQVKIIVPNSTAGLIIGKGGATVKAVMEQSGAWVQLSQKPDGINLQERVTVSGEPEQNRKAV ELIIQKIQEDPQSGSCLNISYANVTGPVANSNPTGSPYANTAEVLPTAAAAAGLLGHANLAGVAAFPVAVL SGFTGNDLVAITSALNLTASYGYNLNTLGLGLSQAATGALAAAAASANPAAAAANLLATYASEASASGS TAGGTAGTFALGSLAAATAATNGYFGAASPLAASAILGTEKSTDGSKDVVEIAPENLVGAILGKGGKTL VEYQELTGARIQISKKGFEVPGTRNRKVTITGTPAATQAAQYLITQRITYEQGVRAANPQKVG</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	49.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: [NP_006480](#)

Locus ID: 4857

UniProt ID: [P51513](#)

RefSeq Size: 3846


Cytogenetics: 14q12

RefSeq ORF: 1449

Synonyms: Nova-1

Summary: This gene encodes a neuron-specific RNA-binding protein, a member of the Nova family of paraneoplastic disease antigens, that is recognized and inhibited by paraneoplastic antibodies. These antibodies are found in the sera of patients with paraneoplastic opsoclonus-ataxia, breast cancer, and small cell lung cancer. Alternatively spliced transcripts encoding distinct isoforms have been described. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified NOVA1 protein (Cat# TP321328). The protein was produced from HEK293T cells transfected with NOVA1 cDNA clone (Cat# [RC221328]) using MegaTran 2.0 (Cat# [TT210002]).