

Product datasheet for **TP321343L**

NXN (NM_022463) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human nucleoredoxin (NXN), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221343 representing NM_022463 Red =Cloning site Green =Tags(s)

MSGFLEELLGEKLVTTGGGEEVDVHSLGARGISLLGLYFGCSLSAPCAQLSASLAIFYGRLRGDAAAGPGP
GAGAGAAAEPEPRRRLEIVFVSSDQDQRQWQDFVRDMPWLALPYKEKHKRLKLWNKYRISNIPSLIFLDA
TTGKVVCNRNGLLVIRDDPEGLEFPWGPKEFREVIAGPLLRNNGQSLESSLEGSHVGVYFSAHWCPPCRS
LTRVLVESYRKIKEAGQNFEIIFVSADRSEESFKQYFSEMPWLAVPYTDEARRSLNRLYGIQGIPTLIM
LDPQGEVITRQGRVEVLNDEDCREFPWHPKPVLELSDSNAAQLNEGPCLVLFVDSSEDDGESEAAKQLIQP
IAEKIIAKYKAKEEEAPLLFFVAGEDDMTDSLRYDTNLPEAAPLLTILDM SARAKYVMDVEEITPAIVEA
FVNDFLAEKLKPEPI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	48.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_071908</u>



[View online »](#)

Locus ID: 64359

UniProt ID: [Q6DKJ4](#)

RefSeq Size: 3004

Cytogenetics: 17p13.3


RefSeq ORF: 1305

Synonyms: NRX; RRS2; TRG-4

Summary: This gene encodes a member of the thioredoxin superfamily, a group of small, multifunctional redox-active proteins. Members of this family are characterized by a conserved active motif called the thioredoxin fold that catalyzes disulfide bond formation and isomerization. The encoded protein acts a redox-dependent regulator of the Wnt signaling pathway and is involved in cell growth and differentiation. [provided by RefSeq, Sep 2015]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified NXN protein (Cat# [TP321343]). The protein was produced from HEK293T cells transfected with NXN cDNA clone (Cat# [RC221343]) using MegaTran 2.0 (Cat# [TT210002]).