

Product datasheet for TP321385L

IREB2 (NM_004136) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human iron-responsive element binding protein 2 (IREB2), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221385 representing NM_004136 Red=Cloning site Green=Tags(s)

MDAPKAGYAFEYLIETLNDSSHKFFDVSKLGTKYDVLPSYRVLLEAAVRNCDGFLMKKEDVMNILDWK
TKQSNVEVPPFFPARVLLQDFTGIPAMVDFAAMREAVKTLGGDPEKVHPACPTDLTVDHSLQIDFSKCAIQ
NAPNPGGGDLQKAGKLSPLKVQPKKLPCRQTTTCRGS CD SGELGRNSGTFSSQIENTPILCFHLQPVPE
PETVLKNQEVEFGRNRERLQFFKWSSRVFKNAVIPPPTGMAHQINLEYLSRVFEEKDLLFPDSVVGTD
SHITMVNGLGILGWGVGGIETEAVMLGLPVSLTLPEVVGCELTGSSNPFVTSIDVWLGITKHLRQVGVAG
KFVEFFGSGVSQLSIVDRTTIANMCPEYGAILSFFPVDNVTLKHLEHTGFSKAKLESMETYLKAVKLFRN
DQNSSGEPEYSQVIQINLNSIVPSVSGPKRPQDRVAVTDMKSDFQAACNEKVGFGKFQIAAEKQKDIVSI
HYEGSEYKLSHGVSVIAAVISCTNNCNPSVMLAAGLLAKKAVEAGLRVKPYIRTSLSPGSGMVTHYLSSS
GVLPLYSLKLGFEIVGYGCSTCVGNTAPLSDAVLNAVKQGDVLTGILSGNKNFEGRLCDCVRANYLASPP
LVWAYAIAGTVNIDFQTEPLGTDPTGKNIYLHDIWPSREEVHRVEEHHVILSMFKALKDKIEMGNKRWNS
LEAPDSVLFPWDLKSTYIRCPFFDKLTKEPIALQAIENAHVLLYLGDSVTTDHPAGSIARNSAAAKY
LTNRGLTPREFNSYGARRGNDVMTGRGTFANIKLFNKFVGPAPKTIHFPSGQTLDFVFAAEALYQKEGIP
LIILAGKKYGSGNSRDWAAKGPYLLGVKAVLAESYEKIHKDHILIGIAPLQFLPGENADSLGLSGRETF
SLTFPEELSPGITLNIQTSTGKVFSVIASFEDDVEITLYKHGGLLN FVARKFS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	104.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004127
Locus ID:	3658
UniProt ID:	P48200 , D3DW85
RefSeq Size:	3928
Cytogenetics:	15q25.1
RefSeq ORF:	2889
Synonyms:	ACO3; IRE-BP 2; IRE-BP2; IRP2; IRP2AD; NDCAMA
Summary:	The protein encoded by this gene is an RNA-binding protein that acts to regulate iron levels in the cells by regulating the translation and stability of mRNAs that affect iron homeostasis under conditions when iron is depleted. When iron levels are low, this protein binds to iron-responsive elements (IRES), stem-loop structures located either in the 5' or 3' UTRs. Binding to the 5' UTR represses translation, while binding to the 3' UTR inhibits mRNA degradation. When iron is found in the cell, this protein is degraded in a F-box and leucine rich repeat protein 5-dependent manner. Variants in this gene have been associated with lung cancer and chronic obstructive pulmonary disease (COPD). Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2017]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified IREB2 protein (Cat# [TP321385]). The protein was produced from HEK293T cells transfected with IREB2 cDNA clone (Cat# [RC221385]) using MegaTran 2.0 (Cat# [TT210002]).