

Product datasheet for TP321398L

FAK (PTK2) (NM_005607) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Recombinant protein of human PTK2 protein tyrosine kinase 2 (PTK2), transcript variant 2, 1 mg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC221398 representing NM_005607
Red=Cloning site Green=Tags(s)

MISADCNLCLEPYDRYLASSKIMAAAYLDPNLNHTPNSSKTHLGTGMERSPGAMERVLKVFHYFESNSE
 PTTWASIIIRHGDATDVRGIIQKIVDSHKVVKHVACYGFRLSHLRSEEVHVLHVDMGVSSVREKYELAHPPE
 EWKYELRIRYLPKGFLNQFTEDKPTLNFFYQVQKSDYMLEIADQVDQEIALKLGCLEIRRSYWEMRGNAL
 EKKSNEYEVLEKDVGLKRFFPKSLDSVKAKTLRKLIIQQTFRQFANLNREESILKFFEILSPVYRFDKECF
 KCALGSSWIIISVELAIGPEEGISYLTDKGCNPTHLADFTQVQTIQYSNSEDKDRKGMLQLKIAGAPEPLT
 VTAPSLTIAENMADLIDGYCRLVNGTSQSFIIRPQKEGERALPSIPKLANSEKQGMRTHAVSVSETDDYA
 EIIDEEDTYTMPSTRDYEIQRERIELGRCIGEGQFGDVHQGIYMSPENPALAVAIKTCCKNCTSDSVREKF
 LQEALTMRFQDHPHIVKLVITENPVWIIMELCTLGELRSFLQVRKYSLDLASELILYAYQLSTALAYLE
 SKRFVHRDIAARNVLSNDCVKLGDFGLSRYMEDSTYYKASKGKLPKWMAPESINFRFRFTSASDVWWMF
 GVCMWELMHGVPFQGVKNNDVIGRIENGERLPMPPNCPPTLYSLMTCWAYDPSRRRPFTELKAQLST
 ILEEEKAQQEERMESRRQATVSWDSGGSDAPPKPSRPGYSPRSSEGFYSPQHMVQTNHYQVSGYP
 GSHGITAMAGSIYPGQASLLDQTDSWNHRPQEIAMWQPNVEDSTVLDLRGIGQVLPHTLMEERLIRQQQE
 MEEDQRWLEKEERFLKPDVRLSRGSDREDGSLQGPIGNQHIYQPVGKPDPAAPPKPPRPGAPGHLGSL
 ASLSSPADSYNEGKLVKQPQEIAPPPTANLDRSNDKVYENVTLVKAVIEMSSKIQPAPPEEYVPMVKEVG
 LALRTLLATVDETIPLLPASTHREIEMAQKLLNSDLGELINKMKLAQQYVMTSLQQEYKKQMLTAAHALA
 VDAKNLLDVIDQARLKMLGQTRPH


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 121.5 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_005598
Locus ID:	5747
UniProt ID:	Q658W2 , Q59GM6
RefSeq Size:	4442
Cytogenetics:	8q24.3
RefSeq ORF:	3222
Synonyms:	FADK; FADK 1; FAK; FAK1; FRNK; p125FAK; pp125FAK; PPP1R71
Summary:	This gene encodes a cytoplasmic protein tyrosine kinase which is found concentrated in the focal adhesions that form between cells growing in the presence of extracellular matrix constituents. The encoded protein is a member of the FAK subfamily of protein tyrosine kinases but lacks significant sequence similarity to kinases from other subfamilies. Activation of this gene may be an important early step in cell growth and intracellular signal transduction pathways triggered in response to certain neural peptides or to cell interactions with the extracellular matrix. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2017]
Protein Families:	Druggable Genome, Protein Kinase
Protein Pathways:	Axon guidance, Chemokine signaling pathway, ErbB signaling pathway, Focal adhesion, Leukocyte transendothelial migration, Pathways in cancer, Regulation of actin cytoskeleton, Small cell lung cancer, VEGF signaling pathway

Product images:

Coomassie blue staining of purified PTK2 protein (Cat# [TP321398]). The protein was produced from HEK293T cells transfected with PTK2 cDNA clone (Cat# [RC221398]) using MegaTran 2.0 (Cat# [TT210002]).