

## Product datasheet for **TP321441L**

### **C15ORF27 (TMEM266) (NM\_152335) Human Recombinant Protein**

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 15 open reading frame 27 (C15orf27), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221441 representing NM_152335 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MAVAPSFNMTNPQPAIEGGISEVEIISQQVDEETKSIAPVQLVNFAYRDLPLAAVDLSTAGSQLLSNLDE DYQREGSNWLKPCCGKRAAVWQVFLLSASLNSFLVACVILWILLTLELLIDIKLLQFSSAFQFAGVIHW ISLVILSVFFSETVLRIVVLGIWDYIENKIEVFDGAVIILSLAPMVASTVANGPRSPWDAISLIIMLRIW RVKRVIDAYVLPVKLEMEMVIQYQYKAKVIQDEQLERLTQICQEQQFEIRQLRAHLAQDDLAAEREA LQAPHVLSQPRSRFKVLEAGTWDEETAESVVEELQPSQEATMKDDMNSYISQYYNGPSSDSGVPEPAVC MVTAAIDIHQPNISSDLFSLDMPLKLGNGTSATSEASRSSHVTRAQSDSSQTLGSSMDCSTAREEPSS EPGSPPLPSQQVVEATVQDLLSLSLSEDPSPQKALDPAPLARSPAGSAQTSPELEHRVSLFNQKNQ EGFTVFQIRPVIHFQPTVPMLEDKFRSLESKEQLHRVPEA  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	58.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online >](#)

RefSeq: [NP\\_689548](#)

Locus ID: 123591

UniProt ID: [Q2M3C6](#)

RefSeq Size: 2414

Cytogenetics: 15q24.2


RefSeq ORF: 1593

Synonyms: C15orf27; HsHVRP1; hTMEM266; HVRP1

**Summary:** Voltage-sensor protein present on the post-synaptic side of glutamatergic mossy fibers and granule cells in the cerebellum (PubMed:25165868, PubMed:30810529). Despite the presence of a voltage-sensor segment, does not form a functional ion channel and its precise role remains unclear (PubMed:25165868, PubMed:30810529). Undergoes both rapid and slow structural rearrangements in response to changes in voltage (PubMed:30810529). Contains a zinc-binding site that can regulate the slow conformational transition (PubMed:30810529). [UniProtKB/Swiss-Prot Function]

**Protein Families:** Druggable Genome, Transmembrane

### Product images:



Coomassie blue staining of purified TMEM266 protein (Cat# [TP321441]). The protein was produced from HEK293T cells transfected with TMEM266 cDNA clone (Cat# [RC221441]) using MegaTran 2.0 (Cat# [TT210002]).