

Product datasheet for TP321446M

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

alpha Synuclein (SNCA) (NM_007308) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human synuclein, alpha (non A4 component of amyloid precursor)

(SNCA), transcript variant NACP112, 100 μg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC221446 representing NM_007308

or AA Sequence: Red=Cloning site Green=Tags(s)

MDVFMKGLSKAKEGVVAAAEKTKQGVAEAAGKTKEGVLYVGSKTKEGVVHGVATVAEKTKEQVTNVGGAV

VTGVTAVAQKTVEGAGSIAAATGFVKKDQLGKEGYQDYEPEA

TRTRPLEQKLISEEDLAANDILDYKDDDDK**V**

Tag: C-Myc/DDK

Predicted MW: 11.2 kDa

Concentration: >0.05 μg/μL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 009292

Locus ID: 6622

UniProt ID: <u>P37840</u>, <u>H6UYS5</u>

RefSeq Size: 1096



Cytogenetics: 4q22.1

RefSeq ORF: 336

Synonyms: NACP; PARK1; PARK4; PD1

Summary: Alpha-synuclein is a member of the synuclein family, which also includes beta- and gamma-


synuclein. Synucleins are abundantly expressed in the brain and alpha- and beta-synuclein inhibit phospholipase D2 selectively. SNCA may serve to integrate presynaptic signaling and membrane trafficking. Defects in SNCA have been implicated in the pathogenesis of Parkinson disease. SNCA peptides are a major component of amyloid plaques in the brains of patients with Alzheimer's disease. Alternatively spliced transcripts encoding different isoforms have

been identified for this gene. [provided by RefSeq, Feb 2016]

Protein Families: Druggable Genome

Protein Pathways: Alzheimer's disease, Parkinson's disease

Product images:



Coomassie blue staining of purified SNCA protein (Cat# [TP321446]). The protein was produced from HEK293T cells transfected with SNCA cDNA clone (Cat# [RC221446]) using MegaTran 2.0 (Cat# [TT210002]).