

Product datasheet for **TP321495M**

PTPRS (NM_130854) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protein tyrosine phosphatase, receptor type, S (PTPRS), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T



[View online »](#)

Expression cDNA Clone >RC221495 representing NM_130854
or AA Sequence: Red=Cloning site Green=Tags(s)


MAPTWGPGMVSVVGP MGLLVLLVGGCAAEEPPRFIPEPKDQIGVSGGVASFVCQATGDPKPRVTWNKKG
 KKVNSQRFETIEFDESAGAVLRIQLRTPRDENVYECVAQNSVGEITVHAKLTVLREDQLPSGFNPIDMG
 PQLKVVERTRTATMLCAASGNPDPEITWFKDFLPVDPSASNGRIKQLRSGALQIESSEETDQGKYECVAT
 NSAGVRYSSPANLYVRRVAPRFSILPMSHEIMPGGNVITCVAVGSPMPYVKWMQGAEDLTPEDDMPV
 GRNVLELTDVKDSANYTCVAMSSLGVIEAVAQITVKSLPKAPGTPMVTENTATSITITWDSGNPDVSVY
 VIEYKSKSQDGPYQIKEDITTTYSIGGLSPNSEYEIWVSAVNSIGQGPPSESVTRTGEQAPASAPRNV
 QARMLSATTMIVQWEEPVEPNGLIRGYRVYTMPEHPVGNWQKHNVDSDLLTTVGSLLLEDETYTVRVLV
 FTSVGDGPLSDPIQVKTQQGVPQPMNLRAEARSETSITLSWSPRQESIIKYELLFREGDHGREVGRTF
 DPTTSYVWEDLKNTEYAFRLAARSPQGLGAFTPVVRQRTLQSKPSAPPQDVKCVSRSTAILVSWRPPP
 PETHNGALVGYSVRYRPLGSEDPEPKEVNGIPPTTTQILLEALEKWTQYRITVAHTEVGPPESSPVV
 RTDEDVPSAPPRKVEAEALNATAIRVLWRSPAPGRQHGQIRGYQVHYVRMEGAEARGPPRIKDVMLADAQ
 EMVITNLQPETAYSITVAAYTMKGDGARSKPKVWVTKGAVLGRPTLSVQQTPEGSLLARWEPPAGTAEDQ
 VLG YRLQFGREDSTPLATLEFPSED RYTAGSVHKGATYVFRLAARSRGLGEEAAEVLSPEDTPRGHP
 QILEAAGNASAGTVLLRWLPPVPAERNGAIVKYTVAVREAGALGPARETELPAEAEPGAENALTLQGLKP
 DTAYDLQVRAHTRRGPGPFSPVRYRTFLRDQVSPKNFKVKMIMKTSVLLSWEFPDNYNSPTPYKIQYNG
 LTLDVDGRTTKLLITHLKPHTFYNFVLTNRGSSLGGLQQTVAWTA FNLLNGKPSVAPKPADGFMVYL
 PDGQSPVPVQSYFIVMVPLRKSRRGGQFLPLGSPEDMDLEELIQDISRLQRRSLRHSRQLEVPRPYAAR
 FSVLPPTFHPGDQKQYGGFDNRGLEP GHRYVFLVLAFLVQKSEPTFAASPFSDPFQLDNPDPQPIVDGEEG
 LIWVIGPVLAVFIICIVIAILLYKNKPDSKRKDSEPRTKCLLNADLAPHHPKDPVEMRRINFQTPGML
 SHPIPIADMAEHTERL KANDSLKLSQEYESIDPGQQTWEHSNLEVNKPKNRYANVIAYDHSRVLQPI
 EGIMGSDYINANYVDGYRRQNAYIATQGPLPETFGDFWRMVWEQRSATIVMMTRLEEKSRKCDQYWPNR
 GTETYGFIQVTLTDLTIELATFCVRTFSLHKNGSSEKREVRQFQFTAWPDHGVPEYPTPFLAFLRRVKTCN
 PPDAGPIVHCSAGVGRGTCFIVIDAMLERIKPEKTVDVYGHVTLMRSQRNYMVQTEDQYSFIHEALLEA
 VGCNGTEVPARSLYAYIQKLAQVEPGEHVTGMELEFKRLANSKAHTSRFISANLPCNKFKNRLVNIMPYE
 STRVCLQPIRGVEGSDYINASFIDGYRQQKAYIATQGPLAETTEDFWRMLWENNSTIVMLTKLREMGRE
 KCHQYWPAERSARYQYFVVDPMAEYNMPQYILREFKVTDARDGQSRTVRQFQFTDWPEQGVPKSGEGFID
 FIGQVHKTKEQFGQDGPISVHCSAGVGRGTFITLSIVLERMRYEGVVDIFQTVKMLRTQRPAMVQTEDE
 YQFCYQAALYLG SFDHYAT

SGPTRRRRLEQKLISEEDLAANDILDYKDDDDKV

- Tag:** C-Myc/DDK
- Predicted MW:** 209.7 kDa
- Concentration:** >0.05 µg/µL as determined by microplate BCA method
- Purity:** > 80% as determined by SDS-PAGE and Coomassie blue staining
- Buffer:** 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
- Preparation:** Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_570924
Locus ID:	5802
UniProt ID:	Q13332 , Q8NHSZ , Q59FX6
RefSeq Size:	6386
Cytogenetics:	19p13.3
RefSeq ORF:	5730
Synonyms:	PTPSIGMA; R-PTP-S; R-PTP-sigma
Summary:	The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP contains an extracellular region, a single transmembrane segment and two tandem intracytoplasmic catalytic domains, and thus represents a receptor-type PTP. The extracellular region of this protein is composed of multiple Ig-like and fibronectin type III-like domains. Studies of the similar gene in mice suggested that this PTP may be involved in cell-cell interaction, primary axonogenesis, and axon guidance during embryogenesis. This PTP has been also implicated in the molecular control of adult nerve repair. Four alternatively spliced transcript variants, which encode distinct proteins, have been reported. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Phosphatase, Transmembrane

Product images:



Coomassie blue staining of purified PTPRS protein (Cat# [TP321495]). The protein was produced from HEK293T cells transfected with PTPRS cDNA clone (Cat# [RC221495]) using MegaTran 2.0 (Cat# [TT210002]).