

Product datasheet for TP321498

Protein Kinase D2 (PRKD2) (NM_001079881) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protein kinase D2 (PRKD2), transcript variant 3, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221498 representing NM_001079881 Red=Cloning site Green=Tags(s)

MATAPSYAGLPGSPGPGSPPPPGGLELQSPPLLPQIPAPGSGVSFHIQIGLTREFVLLPAASELAHVK
QLACSIVDQKFPECGFYGLYDKILLFKHDPTSANLLQLVRSSGDIQEGDLVEVWLSASATFEDFQIRPHA
LTVHSYRAPAFCDHCGEMLFGLVRQGLKCDGCLNYHKRCAFSIPNNCSGARKRRLSSTSLASGHSVRLG
TSESLPCTAEELSRSTTELLPRRPPSSSSSSASSYTGRPIELDKMLLSKVKVPHTFLIHSYTRPTVCQA
CKKLLKGLFRQGLQCKDCKFNCHKRCATRVPNDCLEALINGDVPMEEATDFSEADKSALMDESEDSGVI
PGSHSENALHASEEEEEGEGGKAQSSLGYIPLMRVVQSVRHTRKSTTLREGWVWHYSNKDTLRKRHYWR
LDCKITLQFNNTTNRYYKEIPLSEILTVESAQNFLVPPGTNPHCFEIVTANATYFVGEMPPGGTTPGGPS
GQGAEAARGWETAIRQALMPVILQDAPSAPGHAPHRQASLSISVSNSQIQENVDIATVYQIFPDEVLGSG
QFGVVYGGKHKRKTGRDVAVKVIDKLRFPKQESQLRNEVAILQSLRHPGIVNLECMFETPEKVFVMEKL
HGDMLEMILSSEKGRLEPERLTKFLITQILVALRHLHFKNIVHCDLKPENVLLASADPPQVKLCDFGFAR
IIGKSFRRSVVGTAPYLAPEVLLNQGYNRSLDMWSVGVIMYVSLSGTFPFNEDEDINDQIQNAAFMYPA
SPWSHISAGAILINLLQVKMRKRYSDKSLSHPWLQEYQTWLDLRELEGKMGERYITHESDDARWEQF
AAEHPLPGSGLPTDRDLGGACPPQDHDMQGLAERISVL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	96.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_001073350](#)

Locus ID: 25865

UniProt ID: [Q9BZL6](#)

RefSeq Size: 3202

Cytogenetics: 19q13.32


RefSeq ORF: 2634

Synonyms: HSPC187; nPKC-D2; PKD2

Summary: The protein encoded by this gene belongs to the protein kinase D (PKD) family of serine/threonine protein kinases. This kinase can be activated by phorbol esters as well as by gastrin via the cholecystokinin B receptor (CCKBR) in gastric cancer cells. It can bind to diacylglycerol (DAG) in the trans-Golgi network (TGN) and may regulate basolateral membrane protein exit from TGN. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Protein Kinase

Product images:



Coomassie blue staining of purified PRKD2 protein (Cat# TP321498). The protein was produced from HEK293T cells transfected with PRKD2 cDNA clone (Cat# [RC221498]) using MegaTran 2.0 (Cat# [TT210002]).