

Product datasheet for TP321572L

Olfactory Marker Protein (OMP) (NM_006189) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human olfactory marker protein (OMP), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221572 representing NM_006189 Red =Cloning site Green =Tags(s) MAEDRPQQPQLDMPLVLDQGLTRQMRLRVESLKQRGEKRQDGEKLLQPAESVYRLNFTQQRLQFERWNV VLDKPGKVTITGTSQNWTPDLTNLMTRQLLDPTAIFWRKEDSDAIDWNEADALEFGERLSDLAKIRKVMY FLVTFGEGVEPANLKASVFNQL SGPTRRRLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	18.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006180
Locus ID:	4975
UniProt ID:	P47874
RefSeq Size:	492



[View online »](#)


Cytogenetics: 11q13.5

RefSeq ORF: 489

Summary: Olfactory marker protein is uniquely associated with the mature olfactory receptor neurons in many vertebrate species from fish to man. The OMP gene structure and protein sequence are highly conserved between mouse, rat and human. Results of the mouse knockout studies show that OMP-null mice are compromised in their ability to respond to odor stimuli, and that OMP represents a novel modulatory component of the odor detection/signal transduction cascade. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified OMP protein (Cat# [TP321572]). The protein was produced from HEK293T cells transfected with OMP cDNA clone (Cat# [RC221572]) using MegaTran 2.0 (Cat# [TT210002]).