

Product datasheet for TP321599

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CD1 (CD1A) (NM_001763) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human CD1a molecule (CD1A), 20 μg

Species: Human
Expression Host: HEK293T

Expression cDNA >RC221599 representing NM_001763
Clone or AA Red=Cloning site Green=Tags(s)

Sequence:

MLFLLLPLLAVLPGDGNADGLKEPLSFHVIWIASFYNHSWKQNLVSGWLSDLQTHTWDSNSSTIVFLWPW SRGNFSNEEWKELETLFRIRTIRSFEGIRRYAHELQFEYPFEIQVTGGCELHSGKVSGSFLQLAYQGSDF VSFQNNSWLPYPVAGNMAKHFCKVLNQNQHENDITHNLLSDTCPRFILGLLDAGKAHLQRQVKPEAWLSH GPSPGPGHLQLVCHVSGFYPKPVWVMWMRGEQEQQGTQRGDILPSADGTWYLRATLEVAAGEAADLSCRV

KHSSLEGQDIVLYWEHHSSVGFIILAVIVPLLLLIGLALWFRKRCFC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 35.2 kDa

Concentration: >0.05 μg/μL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling

conditions. Avoid repeated freeze-thaw cycles.

RefSeq: <u>NP 001754</u>

Locus ID: 909



CD1 (CD1A) (NM_001763) Human Recombinant Protein - TP321599

UniProt ID: P06126 2072 RefSeq Size: Cytogenetics: 1q23.1 RefSeq ORF: 981

CD1; FCB6; HTA1; R4; T6 Synonyms:

Summary: This gene encodes a member of the CD1 family of transmembrane glycoproteins, which are

structurally related to the major histocompatibility complex (MHC) proteins and form


heterodimers with beta-2-microglobulin. The CD1 proteins mediate the presentation of primarily lipid and glycolipid antigens of self or microbial origin to T cells. The human genome contains five CD1 family genes organized in a cluster on chromosome 1. The CD1 family members are thought to differ in their cellular localization and specificity for particular lipid ligands. The protein encoded by this gene localizes to the plasma membrane and to recycling vesicles of the early endocytic system. Alternative splicing results in multiple transcript variants. [provided by

RefSeq, Mar 2016]

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Hematopoietic cell lineage

Product images:



Coomassie blue staining of purified CD1A protein (Cat# TP321599). The protein was produced from HEK293T cells transfected with CD1A cDNA clone (Cat# [RC221599]) using MegaTran 2.0 (Cat# [TT210002]).