

Product datasheet for TP321685M

Semaphorin 5A (SEMA5A) (NM_003966) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human sema domain, seven thrombospondin repeats (type 1 and type 1-like), transmembrane domain (TM) and short cytoplasmic domain, (semaphorin) 5A (SEMA5A), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221685 protein sequence Red =Cloning site Green =Tags(s)
	<p>MKGTCVIAWLFSSLGLWRLAHPEAQGTTQCQRTEHPVISYKEIGPWLREFRAKNAVDFSQTLTFDPGQKEL VVGARNYLFRLQLEDLSLIQAVEWECDEATKKACYSKGSKEECQNYIRVLLVGGDRLFTCGTNAFTPVC TNRSLNLTIEHDQISGMARCPYSPQHNSALLTAGGELYAATAMDFPGRDPAIYRSLGILPPLRTAQYN SKWLNENPNFVSSYDIGNFTYFFRENAVEHDCGKTVFSRAARVCKNDIGGRFLEDTWTTFMKARLNCSR PGEVPFYNELQSTFFLPELDLIYGIFTTNVNSIAASAVCVFNLSAIAQAFSGPFKYQENSRSALWLPYPN PNPHFQCGTVDQGLYVNLTERNLQDAQFILMHEVVQPVTTVPSFMEDNSRFSHVAVDVVQGREALVHII YLATDYGTIKKVRVPLNQTSSSCLLEEIEILFPERRRREPIRSLQILHSQSFLVGLREHVVKIPLKRCQFY RTRSTCIGAQDPYCGWDVVMKKCTSLEESLSMTQWEQSISACPTRNLTVDGHFGVWSPWTPCTHTDGSVA GSCLCRTRSCDSPAPQCGGWQCEGPGMEIANCSRNGGWTPWTSWSPCSTTCGIGFQVRQRSCSNPTPRHG GRVCVGNREERYCNEHLLCPPHMFWTGWGPWERCTAQC GGGIQARRRICENGPDCAGCNVEYQSCNTNP CPELKKTPWTPWTPVNISDNGGHYEQRFRYTCKARLADPNLLEVGRQRIEMRYCSSDGTSGCSTDGLSG DFLRAGRYSAHTVNGAWSAWTSWSQCSRDCSRGIRNRKVCNNPEPKYGGMPCLGPSLEYQECNILPCPV DGVWSCWSPWTKSATCGGGHYMRTRSCSNPAPAYGGDICLGLHTEEALCNTQPCPESWSEWSWSECEA SGVQVRARQCILLFPMGSQCSGNTTESRPCVFDNSFIPEVSVARSSSVEEKRCGEFNMFMIAVGLSSSI LGCLLTLVYTYCQRYQQQSHDATVIHPVSPAPLNTSITNHINKLDKYDSVEAIKAFNKNLILEERNKY FNPHLTGKTYSNAYFTDLNNYDEY</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	118.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining



[View online »](#)

Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_003957</u>
Locus ID:	9037
UniProt ID:	<u>Q13591</u> , <u>X5DR95</u>
RefSeq Size:	11825
Cytogenetics:	5p15.31
RefSeq ORF:	3222
Synonyms:	SEMAF; semF
Summary:	This gene belongs to the semaphorin gene family that encodes membrane proteins containing a semaphorin domain and several thrombospondin type-1 repeats. Members of this family are involved in axonal guidance during neural development. This gene has been implicated as an autism susceptibility gene.[provided by RefSeq, Jan 2010]
Protein Families:	Transmembrane
Protein Pathways:	Axon guidance

Product images:



Coomassie blue staining of purified SEMA5A protein (Cat# [TP321685]). The protein was produced from HEK293T cells transfected with SEMA5A cDNA clone (Cat# [RC221685]) using MegaTran 2.0 (Cat# [TT210002]).