

## Product datasheet for TP321759

### FUSIP1 (SRF10) (NM\_006625) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human FUS interacting protein (serine/arginine-rich) 1 (FUSIP1), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221759 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MSRYLRPPNTSLFVRNVADDTRSEDLRREFGRYGPIVDVYVPLDFYTRRPRGFAYVQFEDVRDAEDALHN LDRKWICGRQIEIQFAQGDRKTPNQMKAKEGRNVYSSRYDDYDRYRRSRRSYERRRSRSDYNYRR SYSPRNSRPTGRPRRSRSHSDNDRPNCSWNTQYSSAYYTSRKI  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	22 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_006616</a>
Locus ID:	10772
UniProt ID:	<a href="#">O75494</a>



[View online »](#)

RefSeq Size:	4076
Cytogenetics:	1p36.11
RefSeq ORF:	549
Synonyms:	FUSIP1; FUSIP2; NSSR; PPP1R149; SFRS13; SFRS13A; SRp38; SRp40; TASR; TASR1; TASR2
Summary:	This gene product is a member of the serine-arginine (SR) family of proteins, which are involved in constitutive and regulated RNA splicing. Members of this family are characterized by N-terminal RNP1 and RNP2 motifs, which are required for binding to RNA, and multiple C-terminal SR/RS repeats, which are important in mediating association with other cellular proteins. This protein interacts with the oncoprotein TLS, and abrogates the influence of TLS on adenovirus E1A pre-mRNA splicing. This gene has pseudogenes on chromosomes 4, 9, 14, 18, and 20. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014]
Protein Families:	Transcription Factors
Protein Pathways:	Spliceosome

### Product images:



Coomassie blue staining of purified SRSF10 protein (Cat# TP321759). The protein was produced from HEK293T cells transfected with SRSF10 cDNA clone (Cat# [RC221759]) using MegaTran 2.0 (Cat# [TT210002]).