

## Product datasheet for TP321771

### CD44 (NM\_000610) Human Recombinant Protein

#### Product data:


<b>Product Type:</b>	Recombinant Proteins
<b>Description:</b>	Recombinant protein of human CD44 molecule (Indian blood group) (CD44), transcript variant 1
<b>Species:</b>	Human
<b>Expression Host:</b>	HEK293T
<b>Expression cDNA Clone or AA Sequence:</b>	>RC221771 representing NM_000610 Red=Cloning site Green=Tags(s)
	<p>MDKFWWHAAWGLCLVPLSLAQIDLNITCRFAGVFHVEKNGRYSISRTEAADLCKAFNSTLPTMAQMEKAL SIGFETCRYGFIEGHVVIPRIHPNSICAANNTGVYILTSNTSQYDTYCFNASAPPEEDCTSVTDLPNAFD GPITITIVNRDGRYVQKGEYRTNPEDIYPSNPTDDDVSSGSSERSSTSGGYIFYTFSTVHPIPEDDSP WITDSTDRIPTLMSTSATATETATKRQETWDFWFLPSESKNHLHTTTQMAGTSSNTISAGWEPNE ENEDERDRHLSFSGSGIDDEDFFISSTISTTPRAFDHTKQNDWTQWNPESHNPVLLQTTTRMTDVRN GTTAYEGNWNPEAHPPLIHHEHHEEEETPHSTSTIQATPSSTTEETATQKEQWFGNRWHEGYRQTPKEDS HSTTGAAAASAHTSHPMQGRTPSPEDSSWTDFFNPISHPMGRGHQAGRRMDMDSSHSITLQPTANPNTG LVEDLDRTGPLSMTTQSNSQSFSTSHGLEEDKDHPTTSTLSSNRNDVTGGRRDPNHSEGSTLLEGY TSHYPHTKESRTFIPVTSAKTGSFGVTAVTVGDSNSNVNRSLSGDQDTFHPSGGSHHTHGESDGHSHGS QEGGANTTSGPIRTPQIPEWLILASLLALILAVNCIANSRRRCGQKKKLVNSGNGAVEDRKPSTGLN GEASKSQEMVHLVKNKESSETPDQFMTAETRNLQNVDKIGV</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
<b>Tag:</b>	C-Myc/DDK
<b>Predicted MW:</b>	79.2 kDa
<b>Concentration:</b>	>50 ug/mL as determined by microplate BCA method
<b>Purity:</b>	> 80% as determined by SDS-PAGE and Coomassie blue staining
<b>Buffer:</b>	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
<b>Bioactivity:</b>	ELISA binding assay (PMID: <a href="#">26140900</a> )
<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Applications:</b>	Cell culture: For culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<u>NP_000601</u>
<b>Locus ID:</b>	960
<b>UniProt ID:</b>	<u>P16070</u>
<b>RefSeq Size:</b>	5748
<b>Cytogenetics:</b>	11p13
<b>RefSeq ORF:</b>	2226
<b>Synonyms:</b>	CDW44; CSPG8; ECMR-III; HCELL; HUTCH-I; IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1
<b>Summary:</b>	The protein encoded by this gene is a cell-surface glycoprotein involved in cell-cell interactions, cell adhesion and migration. It is a receptor for hyaluronic acid (HA) and can also interact with other ligands, such as osteopontin, collagens, and matrix metalloproteinases (MMPs). This protein participates in a wide variety of cellular functions including lymphocyte activation, recirculation and homing, hematopoiesis, and tumor metastasis. Transcripts for this gene undergo complex alternative splicing that results in many functionally distinct isoforms, however, the full length nature of some of these variants has not been determined. Alternative splicing is the basis for the structural and functional diversity of this protein, and may be related to tumor metastasis. [provided by RefSeq, Jul 2008]
<b>Protein Families:</b>	Adult stem cells, Cancer stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Stem cell relevant signaling - DSL/Notch pathway, Transmembrane
<b>Protein Pathways:</b>	ECM-receptor interaction, Hematopoietic cell lineage

### Product images:



Coomassie blue staining of purified CD44 protein (Cat# TP321771). The protein was produced from HEK293T cells transfected with CD44 cDNA clone (Cat# [RC221771]) using MegaTran 2.0 (Cat# [TT210002]).