

Product datasheet for TP321771L

CD44 (NM_000610) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human CD44 molecule (Indian blood group) (CD44), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221771 representing NM_000610 Red=Cloning site Green=Tags(s)
	<p>MDKFWWHAAWGLCLVPLSLAQIDLNITCRFAGVFHVEKNGRYSISRTEAADLCKAFNSTLPTMAQMEKAL SIGFETCRYGFIEGHVVIPRIHPNSICAAANTGVYILTSNTSQYDTYCFNASAPPEEDCTSVTDLPNAFD GPITITIVNRDGRYVQKGEYRTNPEDIYPSNPTDDDVSSGSSERSSTSGGYIFYTFSTVHPIPEDDSP WITDSTDRIPATTLMSTSATATETATKRQETWDWFSWFLPSESKNHLHTTTQMAGTSSNTISAGWEPNE ENEDERDRHLSFSGSGIDDEDFFISSTISTTPRAFDHTKQNQDWTQWNPSHSNPEVLLQTTTRMTDVDRN GTTAYEGNWNPEAHPPLIHHEHHEEEETPHSTSTIQATPSSTTEETATQKEQWFGNRWHEGYRQTPKEDS HSTTGTAASAHTSHPMQGRTPSPEDSSWTDFFNPISHPMGRGHQAGRRMDMDSSHSITLQPTANPNTG LVEDLDRTGPLSMTTQQSNSQSFSTSHGLEEDKDHPTTSTLSSNRNDVTGGRRDPNHSEGSTLLLEGY TSHYPHTKESRTFIPVTSAKTGSFGVTAVTVGDSNSNVNRSLSGDQDTFHPSGGSHHTHGSESDGHSHGS QEGGANTTSGPIRTPQIPEWLILASLLALILAVCIAVNSRRRCGQKKKLVNSGNGAVEDRKPISGLN GEASKSQEMVHLVKNKESSETPDQFMTADETRNLQNVDKIGV</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	79.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	ELISA binding assay (PMID: 26140900)
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000601
Locus ID:	960
UniProt ID:	P16070
RefSeq Size:	5748
Cytogenetics:	11p13
RefSeq ORF:	2226
Synonyms:	CDW44; CSPG8; ECMR-III; HCELL; HUTCH-I; IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1
Summary:	The protein encoded by this gene is a cell-surface glycoprotein involved in cell-cell interactions, cell adhesion and migration. It is a receptor for hyaluronic acid (HA) and can also interact with other ligands, such as osteopontin, collagens, and matrix metalloproteinases (MMPs). This protein participates in a wide variety of cellular functions including lymphocyte activation, recirculation and homing, hematopoiesis, and tumor metastasis. Transcripts for this gene undergo complex alternative splicing that results in many functionally distinct isoforms, however, the full length nature of some of these variants has not been determined. Alternative splicing is the basis for the structural and functional diversity of this protein, and may be related to tumor metastasis. [provided by RefSeq, Jul 2008]
Protein Families:	Adult stem cells, Cancer stem cells, Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS, Stem cell relevant signaling - DSL/Notch pathway, Transmembrane
Protein Pathways:	ECM-receptor interaction, Hematopoietic cell lineage

Product images:



Coomassie blue staining of purified CD44 protein (Cat# [TP321771]). The protein was produced from HEK293T cells transfected with CD44 cDNA clone (Cat# [RC221771]) using MegaTran 2.0 (Cat# [TT210002]).