

Product datasheet for TP321811


C7orf57 (NM_001100159) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 7 open reading frame 57 (C7orf57), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221811 representing NM_001100159 Red =Cloning site Green =Tags(s) MRNTSKELQGATHRYAPCDWYYHVPVKRSEKAVDAPPASQIPGLSNLGDHSENLPGTRRYWIKETDSEY VKLAKQGGRPDLLKHFAPGTRKGSPVAYSLPDWYIHHSKPPTASQQEVRAVSMPDYMVHEEFNPDQAN GS YASRRGPFDFDMKTVWQREAELEKEKKLRLPAIDSKYLSKAGTPLGPKNPAGSRLSFPPVPGQKNSSP TNFSKLISNGYKDEWLQQQRADSDKRTPKTSRASVLSQSPRDLEGPQDAARLQDAEASEGPEDTPESSQ SPEESVSASTPAELK TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	32.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_001093629</u>
Locus ID:	136288


[View online »](#)

UniProt ID: Q8NEG2
 RefSeq Size: 2167
 Cytogenetics: 7p12.3
 RefSeq ORF: 885

Product images:


Coomassie blue staining of purified C7orf57 protein (Cat# TP321811). The protein was produced from HEK293T cells transfected with C7orf57 cDNA clone (Cat# [RC221811]) using MegaTran 2.0 (Cat# [TT210002]).