

Product datasheet for **TP321818L**

Midkine (MDK) (NM_001012334) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens midkine (neurite growth-promoting factor 2) (MDK), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221818 protein sequence Red =Cloning site Green =Tags(s)
	 MQHRGFLLLTLLALLALTSAVAKKKDKVKKGGPGSECAEWAWGPCTPSSKDCGVGFREGTCGAQTQRIRC RVPCNWKKEFGADCKYKFENWGACDGGTGTKVRQGLKKARYNAQCQETIRVTKPCTPKTKAKAKAKKGGK GKD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	15.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001012334
Locus ID:	4192
UniProt ID:	P21741



[View online »](#)

RefSeq Size: 1015

Cytogenetics: 11p11.2


RefSeq ORF: 429

Synonyms: ARAP; MK; NFGF2

Summary: This gene encodes a member of a small family of secreted growth factors that binds heparin and responds to retinoic acid. The encoded protein promotes cell growth, migration, and angiogenesis, in particular during tumorigenesis. This gene has been targeted as a therapeutic for a variety of different disorders. Alternatively spliced transcript variants encoding multiple isoforms have been observed. [provided by RefSeq, Jul 2012]

Protein Families: Druggable Genome, Secreted Protein, Transmembrane

Product images:



Coomassie blue staining of purified MDK protein (Cat# [TP321818]). The protein was produced from HEK293T cells transfected with MDK cDNA clone (Cat# [RC221818]) using MegaTran 2.0 (Cat# [TT210002]).