

## Product datasheet for TP321829L

### SNX5 (NM\_152227) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human sorting nexin 5 (SNX5), transcript variant 1, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221829 protein sequence Red=Cloning site Green=Tags(s)

MAAVPELLQQEEDRSKLRVSVVDLNVDPQLIDIPDALSERDKVKFTVHTKTTLPFQSPFVSVTRQHE  
DFVWLHDTLIETTDYAGLIIPPAPTKPDFDGPREKMQKLGEGEGSMTKEEFAKMKQELEAEYLAVFKKTV  
SSHEVFLQRLSSHPVLSKDRNFHFVLEFDQDLSVRRKNTKEMFGGFFKSVKSADEVLFVGVKVVDDFFE  
QEKNFLLINYYNRIKDSCVKADKMTRSHKNVADDYIHTAACLHSLALEEPTVIKKYLLKVAELFEKLRKVE  
GRVSSDEDLKLTLLRYMLNIEAAKDLLYRRRTKALIDYENSNKALDKARLKSQDVKLAEAHQCECCQKF  
EQLSESAKEELINFKRKRVAEFRKNLIEMSELEIKHARNNVSLQSCIDLFKNN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	46.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_689413</a>
Locus ID:	27131



[View online »](#)

UniProt ID: [Q9Y5X3](#)


RefSeq Size: 2308

Cytogenetics: 20p11.23

RefSeq ORF: 1212

**Summary:** This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein functions in endosomal sorting, the phosphoinositide-signaling pathway, and macropinocytosis. This gene may play a role in the tumorigenesis of papillary thyroid carcinoma. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Sep 2013]

### Product images:



Coomassie blue staining of purified SNX5 protein (Cat# [TP321829]). The protein was produced from HEK293T cells transfected with SNX5 cDNA clone (Cat# [RC221829]) using MegaTran 2.0 (Cat# [TT210002]).