

Product datasheet for TP321938L

PTPRS (NM_130855) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human protein tyrosine phosphatase, receptor type, S (PTPRS), transcript variant 4, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221938 representing NM_130855 Red=Cloning site Green=Tags(s)


MAPTWGPGMVSVVGP MGLLWLLVGGCAAEEPPRFI KEPKDQIGVSGGVASFVCQATGDPKPRVTWNKKG
KKVNSQRFETIEFDESAGAVLRIQPLRTPRDENVYECVAQNSVGEITVHAKLTVLREDQLPSGFNPIDMG
PQLKVVERTRTATMLCAASGNPDPEITWFKDFLPVDPSASNGRIKQLRSGALQIESSEETDQGYECVAT
NSAGVRYSSPANLYVRELREVRVAPRFSILPMSHEIMPGGNVNITCVAVGSPMPYVKWMQGAEDLTPED
DMPVGRNVLELTDVKDSANYTCVAMSSLGVIEAVAQITVKS LPKAPGTPMVTENTATSITITWDSGNPDP
VSYVIEYKSKSQDGPYQIKEDITTTTRY SIGGLSPNSEYEIWWSAVNSIGQGPPSESVVTRTGEQAPASA
PRNVQARMLSATTMIVQWEEPVEPNGLIRGYRVYTTMEPEHPVGNWQKHNVDDSLTTVGSLLEDETYTV
RVLAFSTVGDGPLSDPIQVKTQQGVP GQPMNLRAEARSETSITLSWSPPRQESIIKYELLFREGDHGREV
GRTFDPTTSYVWEDLKPNT EYAFRLAARSPQGLGAFTPVVRQRTLQ SISPKNFVKMIMKTSVLLSWEFP
DNYNSPTPYKIQYNGLTLDVDGR TTKKLITHLKPHTFYNFVLNTRGSSLGGLQQTVAWTA FNLLNGKPS
VAPKPDADGFIMVYLPDGGQSPVPVQSYFIVMVPLRKS RGGQFLTPLGSPEDMDLEELIQDISRLQRRLR
HSRQLEVPRPYIAARFVLPPTFHPGDQKQYGGFDNRGLEPGHRYVLFVLAVLQKSEPTFAASPFSDPFQ
LDNPDPQPIVDGEEGLI WVIGPVLAVVFIICIVIAILLYKNKPDSKRKDSEPRTKLLNNADLAPHHPKD
PVEMRRINFQTPGMLSHPP IADMAEHTERL KANDSLKLSQEYESIDPGQQTWEHSNLEVNKPKNRYA
NVIAVDHSRVLQPIEGIMGSDYINANYVDGYRCQ NAYIATQGPLPETFGDFWRMVWEQRSATIVMMTRL
EEKSRIKCDQYWPNRGTETYGFIQVTL LDTIELATFCVRTFSLHKN GSSEKREVRQFQFTAWPDHGVPEY
PTPFLAFLRRVKT CNPPDAGPIVWHCSAGVGR TGC FIVIDAMLERIKPEKTVDVYGHVTLMRSQRNYMVQ
TEDQYSFIHEALLEAVGCGNTEVPARSLYAIQKLAQVEPGEHVTGMELEFKRLANSKAHTSRFISANLP
CNKFKNRLVNIMPYESTRVCLQPIRGVEGSDYIN ASFIDGYRQQKAYIATQG PLAETTEDFWRMLWENNS
TIVVMLTKLREMGREKCHQYWPAERSARYQYFVVD PMAEYNMPQYILREFKVT DARDGQSRTVRQFQFTD
WPEQGVPKSGEGFIDFIGVHKTKEQFGQDGPISVHCSAGVGR TGVFITLSIVLERMRYEGVVDIFQTVK
MLRTQRPAMVQTEDEYQFCYQA ALEYLGSFDHYAT

SGPTRRRRLEQKLISEEDLAANDILDYKDDDDKV



[View online >](#)

Tag:	C-Myc/DDK
Predicted MW:	166.1 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_570925
Locus ID:	5802
UniProt ID:	Q13332 , Q8NHS7 , Q59FX6 , A0A0A0MR60
RefSeq Size:	5171
Cytogenetics:	19p13.3
RefSeq ORF:	4515
Synonyms:	PTPSIGMA; R-PTP-S; R-PTP-sigma
Summary:	<p>The protein encoded by this gene is a member of the protein tyrosine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including cell growth, differentiation, mitotic cycle, and oncogenic transformation. This PTP contains an extracellular region, a single transmembrane segment and two tandem intracytoplasmic catalytic domains, and thus represents a receptor-type PTP. The extracellular region of this protein is composed of multiple Ig-like and fibronectin type III-like domains. Studies of the similar gene in mice suggested that this PTP may be involved in cell-cell interaction, primary axonogenesis, and axon guidance during embryogenesis. This PTP has been also implicated in the molecular control of adult nerve repair. Four alternatively spliced transcript variants, which encode distinct proteins, have been reported. [provided by RefSeq, Jul 2008]</p>
Protein Families:	Druggable Genome, Phosphatase, Transmembrane

Product images:

Coomassie blue staining of purified PTPRS protein (Cat# [TP321938]). The protein was produced from HEK293T cells transfected with PTPRS cDNA clone (Cat# [RC221938]) using MegaTran 2.0 (Cat# [TT210002]).