

## Product datasheet for TP321978L

### C8orf31 (NM\_173687) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 8 open reading frame 31 (C8orf31), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC221978 representing NM_173687 <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MAEKPHQNSCNSVRQLFKTKQLVTHRDRGSCETHRAQGLLAARTTALQRSPLQQEIWESTTALNLPALAP QGLTAKDAHFLGDTDPIQEGARDHAAGGPFQDRQASVAAQTLWERGQGFSRHHGNHLLYSH
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	14.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_775958</a>
Locus ID:	286122
RefSeq Size:	1876
Cytogenetics:	8q24.3
RefSeq ORF:	396



[View online »](#)

Synonyms: FLJ37131; MGC90155

### Product images:



Coomassie blue staining of purified C8orf31 protein (Cat# [TP321978]). The protein was produced from HEK293T cells transfected with C8orf31 cDNA clone (Cat# [RC221978]) using MegaTran 2.0 (Cat# [TT210002]).