

Product datasheet for **TP322003M**

TTC39C (NM_153211) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tetratricopeptide repeat domain 39C (TTC39C), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222003 protein sequence Red =Cloning site Green =Tags(s)
	<p>MSFGASFVSFLNAMMTFEEEEKMQLACDDLKTEKLCESSEEAGVIETIKNKIKKNVDVRKSAPSMVDRLQR QIIADQCQVYLAVLSFVKQELSAYIKGGWILRKAWKIYNKCYLDINALQELYQKKLTEESLTSDAANDNH IVAEGVSEESLNRKLGAVSFGYGLFHLCSMVPPNLLKIINLLGFPDRLQGLSSLMYASESKDMKAPLA TLALLWYHTVVRPFALDGSNDKAGLDEAKEILLKKEAAYPNSSLFMFFKGRIQRLECQINSALTSFHTA LELAVDQREIQHVCLYEIGWCSMIELNFKDAFDSFERLKNESRWSQCYYAYLTAVCQGATGDVDGAQIVF KEVQKLFKRKNNQIEQFSVKAERFRKQTPTKALCVLASIEVLYLWKALPNCSPNLQRMSSQACHEVDDS SVVGLKYLKLLGAIHKCLGNSEDAVQYFQRAVKDELCRQNNLYVQPYACYELGCLLLDKPETVGRGRALLL QAKEDFSGYDFENRLHVRHAALASLRELVPQ</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	59 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_694943
Locus ID:	125488
UniProt ID:	Q8N584
RefSeq Size:	4900
Cytogenetics:	18q11.2
RefSeq ORF:	1566
Synonyms:	C18orf17; HsT2697

Product images:

Coomassie blue staining of purified TTC39C protein (Cat# [TP322003]). The protein was produced from HEK293T cells transfected with TTC39C cDNA clone (Cat# [RC222003]) using MegaTran 2.0 (Cat# [TT210002]).