

Product datasheet for TP322027

SP3 (NM_003111) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human Sp3 transcription factor (SP3), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222027 representing NM_003111 Red=Cloning site Green=Tags(s)

MTAPEKPKVQEEMAALDVDSGGGGGGGGHGEYLQQQQQHNGAVAAAAAAQDTQPSPLALLAATCSKIG
PPSPGDDEEEAAAAGAPAAAGATGDLASALGGAPNRWEVLSATPTTIKDEAGNLVQIPSAATSSGQYV
LPLQNLQNLQIFSVAPGSDSSNGAVSSVQYQVIPQIQSADGQQVQIGFTGSSDNGGINQESSQIQIIPGS
NQTLLASGTPSANIQNLIPQTGQVQVQVAIGGSSFPQGTQVVANVPLGLPGNITFVPINSVDLDSLGLS
GSSQTMTAGINADGHLINTGQAMDSSDINSERTGERVSPDINETNTDLDLFPVPTSSSSQLPVTIDSTGILQ
QNTNSLTSSGQVHSSDLQGNVIQSPVSEETAQNIQVSTAQPVVQHLQLQESQQPTSQAQIVQGITPQT
IHGVQASGQNISQQALQNLQLQNLNPGTFLIQAQTVTPSGQVQVQVQNLQNLQIQNTAAQQITLT
PVQTLTLGQVAAGGAFTSTPVSLSTGQLPNLQTVTVNSIDSAGIQLHPGENADSPADIRIKEEEDPPEEW
QLSGDSTLNTNDLTHLRVQVDEEGDQHQEGKRLRRVACTCPNCKEGGGRGTNLGKKKQHICHIPGCGK
VYGKTSHLRAHLRWHSGERPFVCNWMYCGKRFRTRDELQRHRRHTHTGEKKFVCEPKSRFMRSDHLAKHI
KTHQNKKGIIHSSSTVLASVEAARDTLITAGGTTLILANIQQGSVSGIGTVNTSATSNDILNTNTEIPLQ
LVTVSGNETME

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	81.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_003102
Locus ID:	6670
UniProt ID:	Q02447 , Q86TP0
RefSeq Size:	3920
Cytogenetics:	2q31.1
RefSeq ORF:	2343
Synonyms:	SPR2
Summary:	This gene belongs to a family of Sp1 related genes that encode transcription factors that regulate transcription by binding to consensus GC- and GT-box regulatory elements in target genes. This protein contains a zinc finger DNA-binding domain and several transactivation domains, and has been reported to function as a bifunctional transcription factor that either stimulates or represses the transcription of numerous genes. Transcript variants encoding different isoforms have been described for this gene, and one has been reported to initiate translation from a non-AUG (AUA) start codon. Additional isoforms, resulting from the use of alternate downstream translation initiation sites, have also been noted. A related pseudogene has been identified on chromosome 13. [provided by RefSeq, Feb 2010]
Protein Families:	Druggable Genome, Transcription Factors

Product images:



Coomassie blue staining of purified SP3 protein (Cat# TP322027). The protein was produced from HEK293T cells transfected with SP3 cDNA clone (Cat# [RC222027]) using MegaTran 2.0 (Cat# [TT210002]).