

Product datasheet for TP322058

SLC34A3 (NM_080877) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens solute carrier family 34 (sodium phosphate), member 3 (SLC34A3)
Species:	Human
Expression Host:	HEK293T
Tag:	C-Myc/DDK
Predicted MW:	63.3 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_543153
Locus ID:	142680
RefSeq Size:	2124
Cytogenetics:	9q34.3
RefSeq ORF:	1797
Synonyms:	HHRH; NPTIIc
Summary:	This gene encodes a member of SLC34A transporter family of proteins, and is expressed primarily in the kidney. It is involved in transporting phosphate into cells via sodium cotransport in the renal brush border membrane, and contributes to the maintenance of inorganic phosphate concentration in the kidney. Mutations in this gene are associated with hereditary hypophosphatemic rickets with hypercalciuria. Alternatively spliced transcript variants varying in the 5' UTR have been found for this gene.[provided by RefSeq, Apr 2010]



[View online »](#)

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified SLC34A3 protein (Cat# TP322058). The protein was produced from HEK293T cells transfected with SLC34A3 cDNA clone (Cat# [RC222058]) using MegaTran 2.0 (Cat# [TT210002]).