

## Product datasheet for TP322092L

### C9orf114 (SPOUT1) (NM\_016390) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 9 open reading frame 114 (C9orf114), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222092 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MAERGRKRPCGPGEHGQRIEWRKWKQKKEEKKKWKDLKLMKKLERQRAQEEQAKRLEEEEEAAEKEDRG            RPYTLSVALPGSILDNAQSPELRTYLAGQIARACAIFCVDEIVWFDEEGQDAKTVEGEFRGVGKKGQACV            QLARILQYLECPQYLKKAFFPKHQDLQFAGLLNPLDSPHHMRQDEESEFREGIVDRPTRPGHGSFVNCG            MKKEVKIDKNLEPGLRVTVRLNQQQHPDCKTYHGKVVSSQDPRTKAGLYWGTYVRLASCLSAVFAEAPFQ            DGYDLTIGTSEKSDVASAQLPNFRHALVVFGLQGLEAGADADPNLEVAEPSVLFDLVNTCPGQGSRT            IRTEEAILISLAALQPGLIQAGARHT</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	41.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_057474</a>
Locus ID:	51490



[View online »](#)

UniProt ID:	<a href="#">Q5T280</a>
RefSeq Size:	4302
Cytogenetics:	9q34.11
RefSeq ORF:	1128
Synonyms:	C9orf114; CENP-32; CENP32; HSPC109
Summary:	Required for association of the centrosomes with the poles of the bipolar mitotic spindle during metaphase (PubMed:20813266, PubMed:25657325). Also involved in chromosome alignment (PubMed:20813266). May promote centrosome maturation probably by recruiting A-kinase anchor protein AKAP9 to centrosomes in early mitosis (PubMed:25657325). Binds specifically to miRNA MIR145 hairpin, regulates MIR145 expression at a postranscriptional level (PubMed:28431233).[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome

### Product images:



Coomassie blue staining of purified SPOUT1 protein (Cat# [TP322092]). The protein was produced from HEK293T cells transfected with SPOUT1 cDNA clone (Cat# [RC222092]) using MegaTran 2.0 (Cat# [TT210002]).