

Product datasheet for **TP322104L**

CCDC32 (NM_052849) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 15 open reading frame 57 (C15orf57), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222104 representing NM_052849 Red =Cloning site Green =Tags(s)
	 MKMFESADSTATRSGQDLWAEICSLPNPEQEDGANNAFSDSFVDSCEGEGQREVADFAVQPAVKPWAP LQDSEVYLASLEKKLRRRIKGLNQEVTSKDMLRTLAQAKKECWDRFLQEKLASEFFVDGLDSESTLEHFK RWLQPKVAVSTEEVQYLIPPESQVEKPVAEDEPAAGDKPAAAEQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	20.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_443081</u>
Locus ID:	90416
UniProt ID:	<u>Q9BV29</u> , <u>A0A024R9J9</u>



[View online »](#)

RefSeq Size: 1724
Cytogenetics: 15q15.1
RefSeq ORF: 555
Synonyms: C15orf57

Product images:

Coomassie blue staining of purified CCDC32 protein (Cat# [TP322104]). The protein was produced from HEK293T cells transfected with CCDC32 cDNA clone (Cat# [RC222104]) using MegaTran 2.0 (Cat# [TT210002]).