

## Product datasheet for TP322159

### URG4 (URGCP) (NM\_001077663) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens up-regulated gene 4 (URG4), nuclear gene encoding mitochondrial protein, transcript variant 3, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222159 representing NM_001077663 Red=Cloning site Green=Tags(s)

MASPGIEVELLGKGHSDLGEVAPEIKASERRTAVAIADLEWREMEGDDCEFRYGDGTNEAQDNDFFPTVER  
SRLQEMLSLGLEYQVQKLSLQDSLQISFDSMKNWAPQVPKDLPNFLRKLQALNADARNTTMVLDVLP  
DARPEKESQMEEIYWDPAADLAADIYSFSELPTDTPVNPLDLLCALLSSDSFLQQEIALKMALCQ  
FALPLVLPDSENHYHTFLLWAMRGIVRTWWSQPPRGMGSFREDSSVLSRAPAFAFVRMDVSSNSKSQLLN  
AVLSPGHRQWDCFWHRDLNLGTNAREISDGLVEISWFFPSGRELDIFPEPVAFLNLRGDIGSHWLQFKL  
LTEISSAVFILTDNISKKEYKLLYSMKESTTKYFILSPYRGKRNTNLRFLNKLIPVLKIDHSHVLVKVS  
STDSDFVKRIRAIVGNVLRAPCRRVSVEDMAHAARKLGLKVDCECQKAKDRMERITRIKIKDSDAYR  
RDELRLQGDPWRKAAQVEKEFCQLQWAVDPPEKHRAELRRRLELRMQQNGHDPSSGVQEFISGISSPSL  
SEKQYFLRWMWGLARVAQPRLRQPPETLLLRPKHGGTTDVGEPLWPEPLGVEHFLREMGQFYEAESCL  
VEAGRLPAGQRRFAHFPGLASELLLTGLPLELIDGSTLSMPVRWVTGLLKEHVRLERRSRLVVLSTVGV  
PGTGKSTLLNTMFGLRFATGKSCGPRGAFMQLITVAEGFSQDLGCDHILVIDSGGLIGGALTSAGDRFEL  
EASLATLLMGLSNVTVISLAETKDIPAAILHAFLRLEKTGHMPNYQFVYQNLHDVSVGPVRPRDKRQLLD  
PPGDLSRAAAQMEKQGDGFRALAGLAFCDPEKQHIWHIPGLWHGAPPMMAAVSLAYSEAIFELKRCLLENI  
RNGLSNQNKNIQQLIELVRRL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	104.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol



[View online »](#)

<b>Preparation:</b>	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
<b>Note:</b>	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
<b>Storage:</b>	Store at -80°C.
<b>Stability:</b>	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
<b>RefSeq:</b>	<a href="#">NP_001071131</a>
<b>Locus ID:</b>	55665
<b>UniProt ID:</b>	<a href="#">Q8TCY9</a>
<b>RefSeq Size:</b>	3636
<b>Cytogenetics:</b>	7p13
<b>RefSeq ORF:</b>	2793
<b>Synonyms:</b>	URG4
<b>Summary:</b>	URG4 is upregulated in the presence of hepatitis B virus (HBV)-encoded X antigen (HBxAg) and may contribute to the development of hepatocellular carcinoma by promoting hepatocellular growth and survival (Tufan et al., 2002 [PubMed 12082552]).[supplied by OMIM, Mar 2008]

### Product images:



Coomassie blue staining of purified URGCP protein (Cat# TP322159). The protein was produced from HEK293T cells transfected with URGCP cDNA clone (Cat# [RC222159]) using MegaTran 2.0 (Cat# [TT210002]).