

Product datasheet for TP322174

EXD3 (NM_017820) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens exonuclease 3'-5' domain containing 3 (EXD3), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222174 representing NM_017820 Red=Cloning site Green=Tags(s)
	<p>MDPGDPAGDPAAGERHRMGQDPLLLLQALQTLWSTREKQLREEAWRGFAALDDPLAGLLDMLESCRGQR GEGPSLAAWISHQLQCWLQAQPCPSLAQHSLRLKQLQARAVKVLTESPPSLAAPLASIFQLQDADRSCLL AHVHRLHHEGRFREAATLGTTLKLQSELGVEKMSIPLLLQDKVALVERYVAGFPDLQRRLLVLMDSWCQP GFDIKDVARRYPEVTSLSLEKLSPKALSRQVLRQLERYGVAPALCPNAAIQRLAALRHLCHKRFVEKSL SQENWTDHVQGLVGQSPWLQEQLSLLVSHSDPVTAQAACAMELLLPEERLPAAVAVELRRFRLQGRATEA DSRLEVKDMKDRYYQLPIPRENVHLLASWEDLTRHEGALLQCHQVVGVDVEWTPVFVAGGRPRPSLLQVA VEGHVFLLDVLAALSQPPTGQGAQAFSRLVAQLLSDPSITKLGYGMMVGDQLKLGTSICALAHVEKQILGGM DLLLVHRQMRVASVPAPAVDRARELRGLSLLVQVQLGTALDKTQQLSNWDRRPLYEEQVIYAAADAYCLL EVHQCALCREPARFHLSDELGSRPRHRERPGARKPPGLQKASAPAAPRQVPVAVAVSEGAAPQIPARAF RVCDNMLQGLARSLRCLGVDARMLGNGEDHRRAAEVARQEGRIILTSGQPFHKLRAQVAGRCLSDVDCS LKAQQQAKAVLKHFNVRVTHADIFSRCQACNCDQYLKVS RDMMKQLMWLSSHQEGPRSSGDEATQSQAVQ EPGPAPDAAPEGCTYDRPCRWLQMA DLRAETPDMLADGTRLQLAGVPVGLRTPGLRCFCYCTGCGKVFV DGSHLGRVATHFRDMLESAPSPCEPSPAPSPASSPF</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	96.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_060290
Locus ID:	54932
UniProt ID:	Q8N9H8
RefSeq Size:	2887
Cytogenetics:	9q34.3
RefSeq ORF:	2628
Synonyms:	mut-7; Nbr
Summary:	Possesses 3'-5' exoribonuclease activity. Required for 3'-end trimming of AGO1-bound miRNAs (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified EXD3 protein (Cat# TP322174). The protein was produced from HEK293T cells transfected with EXD3 cDNA clone (Cat# [RC222174]) using MegaTran 2.0 (Cat# [TT210002]).