

Product datasheet for **TP322177L**

FBXO36 (NM_174899) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens F-box protein 36 (FBXO36), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222177 representing NM_174899 Red =Cloning site Green =Tags(s) MASWLPETLFETVGQGGPPPSKDYYQLLVTRSQVIFRWWKISLRSEYRSTKPGEAKETHEDFLENSHLQGGQ TALIFGARILDYVINFCCKGKFDLFLERLSDDLLLTIISYLDLEDIARLCQTSHRFAKLCMSDKLWEQIVQS TCDTITPDVRLAEDTGWRQLFFTNKLQLRQLRKRKQKYGNLREKQP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	21.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_777559
Locus ID:	130888
UniProt ID:	Q8NEA4
RefSeq Size:	924



[View online »](#)

Cytogenetics: 2q36.3


RefSeq ORF: 564

Synonyms: Fbx36

Summary: Members of the F-box protein family, such as FBXO36, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).[supplied by OMIM, Mar 2008]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified FBXO36 protein (Cat# [TP322177]). The protein was produced from HEK293T cells transfected with FBXO36 cDNA clone (Cat# [RC222177]) using MegaTran 2.0 (Cat# [TT210002]).