

Product datasheet for TP322233

AGO3 (NM_024852) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens eukaryotic translation initiation factor 2C, 3 (EIF2C3), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222233 representing NM_024852 Red=Cloning site Green=Tags(s)

MEIGSAGPAGAQPLLMVPRRPGYGTMGKPIKLLANCFQVEIPKIDVYLYEVDIKPDKCPRRVNREVVDSM
VQHFKVTFGDRRPVYDGRSLYTANPLPVATTGVDLDTLPGEKDRPFKVSIFVSRVSWHLLHEVL
TGRTLPEPLELDKPISTNPVHAVDVLRHPLSMKYTPVGRSFFSAPEGYDHPPLGGGREVWFGFHQSVRPA
MWKMMLNIDVSATAFYKAQPVIQFMCEVLDIHNIDEQPRPLTDSHRVKFTKEIKGLKVEVTHCGTMRRKY
RVCNVTRRPASHQTFPLQLENGQTVERTVAQYFREKYTLQKYPHLPCLQVGEQKHTYLPLEVCNIVAG
QRCIKKLTNDQSTMIKATARSAPDRQEEISRLVRSANYETDPFVQEFQFKVRDEMAHVTGRVLPAPMLQ
YGGNRNRTVATPSHGVDWDMRGKQFHTGVEIKMWAIACFATQRQCEEILKGFTDQLRKISKDAGMPIQQQP
CFCKYAQGADSVPEPMFRHLKNTYSGLQLIIVILPGKTPVYAEVKRVGDTLLGMATQCVQVKNVIKTSPQT
LSNLCLKINVKLGGINNILVPHQRPSVFQQPVI FLGADVTHPPAGDGKKPSIAAVVGSMDAHPSTRYCATV
RVQRPRQEIIQDLASMVRELLIQFYKSTRFKPTRIIFYRDGVSEGGFRQVLYYELLAIREACISLEKDYQ
PGITYIVVQKRHHTRLFCADRTERVGRSGNIPAGTTVDTDITHPYEFDYLCSHAGIQGTSRPSHYHVLW
DDNCFTADELQLLTYQLCHTYVRCTRSVSIPAPAYYAHVAFRARYHLVDKEHDSAEGSHVSGQSNRDP
QALAKAVQIHQDTLRMTMYFA

SGPTRRRLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	97.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.



[View online »](#)

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_079128](#)

Locus ID: 192669

UniProt ID: [Q9H9G7](#), [B4E1P5](#)

RefSeq Size: 3578

Cytogenetics: 1p34.3


RefSeq ORF: 2580

Synonyms: EIF2C3

Summary: This gene encodes a member of the Argonaute family of proteins which play a role in RNA interference. The encoded protein is highly basic, contains a PAZ domain and a PIWI domain, and may play a role in short-interfering-RNA-mediated gene silencing. This gene is located on chromosome 1 in a tandem cluster of closely related family members including argonaute 4 and eukaryotic translation initiation factor 2C, 1. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome

Product images:



Coomassie blue staining of purified AGO3 protein (Cat# TP322233). The protein was produced from HEK293T cells transfected with AGO3 cDNA clone (Cat# [RC222233]) using MegaTran 2.0 (Cat# [TT210002]).