

Product datasheet for TP322242

Flt3 ligand (FLT3LG) (NM_001459) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human fms-related tyrosine kinase 3 ligand (FLT3LG), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222242 protein sequence Red =Cloning site Green =Tags(s)

MTVLAPAWSPTTYLLLLLLLSSGLSGTQDCSFQHSPISDFAVKIRELSDYLLQDYPVTVASNLQDEELC
GGLWRLVLAQRWMERLKT VAGSKMQGLLERVNT EIH FVTKCAFQPPPSCLR FVQTNISRLLQETSEQLVA
LKPWITRQNF SRCLELQCQPDSSTLPPWSPRPLEATAPTAPQPPLLLLLLLPVGLLLLAAAWCLHWQRT
RRRTPRPGEQVPPVSPQDLLLLVEH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	26.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Bioactivity:	FLT-3LG generate dendritic cell differentiation: Fms-like tyrosine kinase 3 ligand (Flt3L) is a nonredundant cytokine required for dendritic cell (DC) homeostasis in lymphoid tissues. FLT3L is commonly used to generate pDCs in vitro. In this experiment, 300,000 unfractionated mouse bone marrow cells were cultured with purified recombinant human FLT-3L (OriGene, TP322242) at indicated concentrations. Total live cells were counted at the day 8.
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_001450](#)

Locus ID: 2323

UniProt ID: [P49771](#)

RefSeq Size: 1074

Cytogenetics: 19q13.33

RefSeq ORF: 705


Synonyms: FL; FLG3L; FLT3L


Summary: Dendritic cells (DCs) provide the key link between innate and adaptive immunity by recognizing pathogens and priming pathogen-specific immune responses. FLT3LG controls the development of DCs and is particularly important for plasmacytoid DCs and CD8 (see MIM 186910)-positive classical DCs and their CD103 (ITGAE; MIM 604682)-positive tissue counterparts (summary by Sathaliyawala et al., 2010 [PubMed 20933441]).[supplied by OMIM, Jan 2011]

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein, Transmembrane

Protein Pathways: Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Pathways in cancer

Product images:





Coomassie blue staining of purified FLT3LG protein (Cat# TP322242). The protein was produced from HEK293T cells transfected with FLT3LG cDNA clone (Cat# [RC222242]) using MegaTran 2.0 (Cat# [TT210002]).