

Product datasheet for **TP322278L**

Calpastatin (CAST) (NM_173061) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Purified recombinant protein of Homo sapiens calpastatin (CAST), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC222278 representing NM_173061 Red =Cloning site Green =Tags(s)
	<p>MGQFLSSTFLEGSPATVWHDKLCDGERRGAREAVRIFQDQAKAKEEKLEKCGEDDETIPSEYRLKPATDK DGKPLLPEPEEKPKPRSESELIDELSEDFRSECKEKPSKPTTEESKAAAPAPVSEAVCRTSMCSIQS APPEPATLKGTVPDDAVEALADSLGKKEADPEDGKPVMDKVKEKAKEEDREKLGEKEETIPPDYRLEEVK DKDGKPLLPKESKEQLPPMSEDFLLDALSEDFSGPQNASSLKFEAKLAAISEVVSQTPASTTQAGAPP RDTSSDKDLDDALDKLSDSLGQRQPDENKPMEDKVKEKAKAEHRDKLGERDDTIPPEYRHLLDDNGQD KPKVPPTKKSSEDSKKPADDQDPIDALSGDLDCPSTTETSQNTAKDKCKKAASSSKAPKNGGKAKDSAKT TEETSKPKDD</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	46.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_775084</u>



[View online »](#)

Locus ID:	831
UniProt ID:	P20810
RefSeq Size:	1889
Cytogenetics:	5q15
RefSeq ORF:	1290
Synonyms:	BS-17; calpain inhibitor; calpastatin; heart-type calpastatin; MGC9402; OTTHUMP00000158519; OTTHUMP00000158520; sperm BS-17 component
Summary:	The protein encoded by this gene is an endogenous calpain (calcium-dependent cysteine protease) inhibitor. It consists of an N-terminal domain L and four repetitive calpain-inhibition domains (domains 1-4), and it is involved in the proteolysis of amyloid precursor protein. The calpain/calpastatin system is involved in numerous membrane fusion events, such as neural vesicle exocytosis and platelet and red-cell aggregation. The encoded protein is also thought to affect the expression levels of genes encoding structural or regulatory proteins. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jun 2010]

Product images:



Coomassie blue staining of purified CAST protein (Cat# [TP322278]). The protein was produced from HEK293T cells transfected with CAST cDNA clone (Cat# [RC222278]) using MegaTran 2.0 (Cat# [TT210002]).